NEUATION



DESIGNED TO DELIVER

iStir OH Jr.

(Analog Overhead Stirrer)

NEUATION

iStir OH Jr.

(Analog Overhead Stirrer)

Overhead stirrers are a crucial part of laboratory setups and the work that they do makes the user's life a lot easier. But, in the race to attract customers most manufacturers skimp on the really important stuff and focus on the superficial to sell a lot of devices. In the end, it's the customer who ends up suffering because most of these devices would suffer from some kind of performance or maintenance issues.

Here at Neuation, we start form the basics to create a product that not only performs well, but also creates a pleasant experience for the user. With a complete department focusing on R&D, the designers and product developers at neuation are armed with crucial research about the customer's requirements and needs when it comes to overhead stirrers. The result as you can see is the OH Jr. series of stirrers which combine a compact form factor with exceptional performance for enhanced convenience. Various built-in safety features ensure that no matter what, your safety is never compromised.




NEUATION

Product Features



Brushless DC motor delivers ultra quiet performance and a maintenance-free long life







Variable stirring speed settings upto 2200 RPM for powerful performance

> Built-in overload safety feature to ensure that mindless mistakes do not compromise device & your safety





Small footprint saves valuable bench space

Powerful motor delivers torque upto 10 Ncm to stir solutions with viscosity up to 1000 mPas & 10 L volume



Specifications

iStir OH Jr.		
Specifications	Details	
Motor Type	Brushless DC	
Capacity	10 Liters	
Speed Range	60 - 2200 RPM	
Torque	10 Ncm	
Vascocity	1000 mPas	
Timer	1 min to 99 hr 59 mins & infinity mode	
Protection Class	IP 42	
Dimension (Head)	80 x 42 x 166mm	
Weight	1.3 kgs	

Ordering Information

iStir OH Jr.		
Product Name	Cat #	
iStir OH Jr.	OHS1-1100-EU/US/UK/AU	

Accessories				
Product Name	С	Cat #		
Boss Head	ОН	OHS1-BH		
Impellers	Leng	Length (mm)		
OHS1-IX	300	Jr. impeller		



Neuation_Catalogue_iStir OH Jr._V1