NEUATION



TOP-NOTCH PERFORMER

iFuge M15K

(Table Top Centrifuge)



iFuge M15K

(Table Top Centrifuge)

In any research facility or lab where accuracy reigns supreme, you need an equipment that can deliver precision to the last decimal point. Such an equipment might be hard to get hold of but not impossible, especially when Neutation is here to help you.

With Neuation's iFuge M15K you can expect a lot as it delivers even more. The lightning fast 15000 RPM microprocessor controlled centrifuge can perform at peak efficiency delivering top-notch accuracy in every cycle despite high load.





Witness The Champion



Brushless Dc Motor For Maintanance Free Long Life



Variable Speed From 500-15000 RPM



Microprocessor Controlled Efficient Performance.



Intuitive Interface For Quick And Convenient Setting With Last Run Memory Feature



One Touch RPM / RCF Conversion Operation

Variable Timer Setting From 1 To 99 Min & Infinite Mode



One-touch Short Spin Feature



Easy fit rotor



12 X 1.5 / 2 ml Closed Rotors Design For Lesser Heat Generation



Imbalance Detection With Auto-cut Off Safety Feature



Witness The Champion



Lid-lock Safety Feature Opens The Lid Automatically After Every Run



Enhanced Airflow Design For Less Noise And Greater Sample Safety



Smooth & Silent Operation Even At Maximum Speed Safe Low Voltage Device



Compact Design And Can Also Be Used In Hoods And Cold Rooms



ROTOR OPTION



12 x 1.5/2 ml | 15596 x g



2 x 8 x 0.2 ml strip rotor

Specifications

iFuge M15K		
Specifications	Details	
Rotor Capacity	12 x 1.5ml / 2.0ml Microtubes	
Maximum RPM / RCF	15000 / 15596g	
Display	Large Digital	
Imbalance Detection	Yes	
Motor Type	(BLDC) Brushless DC motor	
Timer Setting	Countdown Timer From 1-99 Mins & Infinity Mode	
Dimensions (WxDxH)	187 x 265 x 115 mm	
Noise Level	< 60 db	
Weight	2.2 Kg	
Rotors and Adaptors (included in standard pack)	Reduction Adaptors For 0.2 ml Microtubes, Reduction Adaptors For 0.5 ml Microtubes	

Ordering Information

iFuge M15K			
Product Name	Cat #	Details	
iFuge M15K	iFuge-M15K UK/AU/US/EU	High Speed Micro Centrifuge Upto 15000 RPM	

Accessories			
Part Name	Cat #	Details	
12 Slot Rotor	M12 R12	12 x 1 5/ 2 ml tubes slot rotor	
PCR Strip Rotor	M12 - R2 x 8	2 PCR stirp rotor	



AN ISO 9001 : 2015 & ISO 13485 : 2016 CERTIFIED