

NEUATION




MORE THAN JUST
ANOTHER ROLLER

iRoll RnR 4
(Four Tube Roller)

iRoll RnR 4 (Four Tube Roller)

Tube rollers are meant to make your work in the lab easier, that's why they come with specifications that advertise high-performance and long durability to justify your investment. But the use of low-quality materials and components to cut production costs results in devices that are a let down in terms of their durability and longevity. Also, due to apparent design flaws, most of these devices are large in size which is a big no-no as the space in your laboratory is already limited and these devices occupy a large space compared to their functionality.

With adequate research into user feedback and their tube rolling requirements, Neuation's R&D department helped integrate valuable insights into the design and development process of iRoll RnR4, resulting in a device which comes with a reliable and durable motor for years of maintenance-free longevity. Keeping in mind the space constraints faced by most users, the exceptional design of iRoll RnR4 manages to wrap incredible performance in a compact package for a wonderful user-experience.



Product Features



Effortless And Tool-free Removal/installation Of Rollers Convenient And Effortless Cleaning.



Easy Adjustable Roller Spacing To Hold Different Size Tubes And Bottles



Fixed Speed Of 30 RPM And Is Accurate Up To +/- 2 RPM.



Gentle Rolling Motion For Optimum Mixing Of Samples With A Tilt Angle Of 3°


Infinite Mode/continuous Mode Of Operation For Enhanced Convenience.



Ac Synchronous Motor Offers



Easy Loading Of Sample While Operational Run



Robust Construction For Stable Operation & Can Be Used In Conditions Upto 40° C And Humidity Upto 80%



Product Features



Compact Design
Makes For A
Small Footprint To Save
Valuable Bench Space.

IP21 Class Compliant
For Added Durability.



Specifications

Specifications	iRoll RnR 4
Capacity	9 x 7 ml tubes
Motor Type	AC Synchronous
Speed	30 RPM
Run Time	Continuous / Infinite
Tilt angle	+/-3°
Speed Accuracy	+/-2 RPM
Noise Level	Noiseless operation
Load Capacity	1.5 kg
Roller Size (mm)	300 (length) x 30 (roller dia) mm
Protection class	IP 21
Dimensions (W X D X H)	475 x 210 x 80 mm
Weight	2 kg

Ordering Information

iRoll RnR 4		
Product Name	Cat #	Details
iRoll RnR 4 (220V)	TR01-3102 EU/UK/AU	4 Tube Roller
iRoll RnR 4 (110V)	TR01-3101-US	4 Tube Roller

Accessories		
Part Name	Cat #	Details
Rollers	TR01 - R1	Extra Roller Attachment
Power Adaptor	TR01 - A1	-

NEUATION

www.neuation.com

AN ISO 9001 : 2015 & ISO 13485 : 2016 CERTIFIED MFG. COMPANY

*We also execute OEM projects. Specifications can be changed without notice for quality improvement.