


# NEUATION



BLEND IT LIKE  
NEVER BEFORE

iSwix FS

Fixed Speed Vortex Mixer

## iSwix FS

### Fixed Speed Vortex Mixer

Vortex mixer are staple mixing systems in laboratory equipment which helps in homogenizing sample of interest.

Much like its Predecessor, with its top of line internal components and class leading feature such as inbuilt counter balance and IP21 rating it delivers superior efficiency enabling accurate and fast homogenization and tackles harder sample matrices easily.



## Product Features



Non Corrosive Body  
For Long And Durable  
Performance

In Built Counter Balance  
Ensures Stable Operation  
For Complete Peace Of Mind



Orbital Diameter  
Of 4 MM

IP21 Compliant Design  
Makes It Spill Safe



AC Synchronous Motor  
For Efficient Noiseless  
Operation

Small Footprint Saves  
Valuable Bench Space And  
Can Also Be Used Inside  
Hoods And Cold Rooms



### Standard Attachments



Rubber Disk  
Attachment



Cup Attachment

### Optional Attachments



Round Disc

# Specifications

iSwix FS	
Max RPM	2800 RPM
Orbital Diameter	4mm
Motor Type	AC Motor
Max Load Capacity	500 gm
Speed Accuracy	± 5%
Mode	Continuous & Touch Mode
Protection Rating	IP 21
Ambient Temperature	5 to 40°C
Dimensions (L X W X H)	205x136x138.5
Weight	3 Kgs
Voltage/Consumption	220V/70W

# Ordering Information

iSwix FS		
Product Name	Cat #	Details
iSwix FS	VM-4210 EU/US/UK/AU	Fixed Speed Vortex Mixer
Attachment	VM-A6	Round Disc
	VM-A7*	Cup Attachment*
	VM-A8*	Rubber Cup Disk Attachment*

\*Supplied as Standard Accessory

## Accessories



Round Disc  
CAT. NO. # VM-A6

# NEUATION

[www.neuation.com](http://www.neuation.com)

AN ISO 9001 : 2015 & ISO 13485 : 2016 CERTIFIED MFG. COMPANY

\*We also execute OEM projects. Specifications can be changed without notice for quality improvement.