

LONG TERM PRESERVATION OF YOUR PRECIOUS SAMPLES


ACCUMAX

Cryovial

(Temperature Stability Upto -196°C/ Volume Range 1/1.8/3.5/4.5ml)



V Bottom For Maximum
Sample Recovery



Autoclavable
121°C For 15 Mins



Better Stability With Hepta
Foot For Self Standing
And Storage Efficiency



Easy Sample Identification,
Large And Clear Patch For
Personalized Labeling




Quick Seal Closure And
Leak Proof Tube Seal



Externally Threaded
Without O-ring To
Avoid Contamination



Certified
DNase RNase Free



For more details visit us at www.accumaximum.com

AN ISO 9001 : 2015 & ISO 13485 : 2016 CERTIFIED MFG. COMPANY

*We also execute OEM projects. Specifications can be changed without notice for quality improvement.