


ACCUMAX



PERFECTION REDEFINED

Accumax Centrifuge Tubes

(15/50 ML Tubes)

## Accumax Centrifuge Tubes


(15/50 ml Tubes)

Where perfection is the prerequisite, even the smallest components get the same importance.

Even though centrifuge tubes are considered widely as just another part of the device, they are extremely important in the way they affect the outcome of the experiment.

Thus, at Accumax, the highest standards are followed during the design & development of these tubes to ensure that there does not remain even the minutest chance of any flaws in this vital consumable.

As a result, these tubes of high quality are certified to be free from all affecting biological elements such as DNase, RNase, Endotoxin, Human DNA and ATP.



## Perfection Is An Understatement

Experience the impact of perfection in your lab with Accumax's impeccable centrifuge tubes which are designed using the absolute best raw materials and procedures to ensure almost 100% purity. Moreover, these tubes are built to offer optimum convenience & durability.



Manufactured In Fully Automatic Human Touch-free Plant With Class 100k Clean Room Facility With Controlled Environment & Access

Transparent And Clear



Produced From Swiss-made Diamond Finish Mold For Flawless 'Super Smooth' Internal Surface

US FDA Approved Virgin Polypropylene Without Any Trace Of Metals



Non Cytotoxic

DNase, RNase, Endotoxin, Human DNA Free



Biologically Inert

Can Withstand Temperature From -80 To +121° C



Trace Metal-Free Tubes Are Also Available

Amber colour tubes ideal for storing light sensitive samples are also available.



Perfection Is An Understatement



Easy Seal Design Of The Cap To Ensure Smooth Opening And Closing With Single Hand Repeatedly



Multiple Packaging Options



Clear Large Graduations On The Body Of The Tube



Superior Engineering Technique Enables Tubes To Withstand Higher G Force

Tamper Proof Graduation Marks & Writing Area - Resistant To Chemical & Physical Corrosion



Fits All Standard And High Capacity Centrifuge Rotors



15 ml  
17000 g

50 ml  
20000 g

# Ordering Information

## Accumax Centrifuge Tubes

Free from Dnase/ Rnase/ Pyrogen/ Endotoxin

Cat. No.	Description	Packing Type	Quantity / Pack	Quantity / Case
ACT-15-B	15 ML Centrifuge Tube With Flat Cap	Bulk	25 Tubes	25 Tubes X 20 packs
ACT-15-R	15 ML Centrifuge Tube With Flat Cap	Rack	25 Tubes	25 Tubes X 20 packs
ACT-15-B-S	15 ML Centrifuge Tube With Flat Cap, Sterile	Bulk	25 Tubes	25 Tubes X 20 packs
ACT-15-R-S	15 ML Centrifuge Tube With Flat Cap, Sterile	Rack	25 Tubes	25 Tubes X 20 packs
ACT-15-S-W	15 ML Centrifuge Tube With Flat Cap, Sterile, Individually Wrapped	Individually Wrapped	25 Tubes	25 Tubes x 20 packs
ACT-15-B-A	15 ML Centrifuge Tube Amber Color With Flat Cap	Bulk	25 Tubes	25 Tubes X 20 packs
ACT-50-B	50 ML Centrifuge Tube With Flat Cap	Bulk	25 Tubes	25 Tubes X 20 packs
ACT-50-R	50 ML Centrifuge Tube With Flat Cap	Rack	25 Tubes	25 Tubes X 20 packs
ACT-50-B-S	50 ML Centrifuge Tube With Flat Cap, Sterile	Bulk	25 Tubes	25 Tubes X 20 packs
ACT-50-R-S	50 ML Centrifuge Tube With Flat Cap, Sterile	Rack	25 Tubes	25 Tubes X 20 packs
ACT-50-S-W	50 ML Centrifuge Tube With Flat Cap, Sterile, Individually Wrapped	Individually Wrapped	25 Tubes	25 Tubes x 20 packs
ACT-50-B-A	50 ML Centrifuge Tube Amber Color With Flat Cap	Bulk	25 Tubes	25 Tubes X 20 packs
ACT-50-K-B-A	50 ml Self Standing Centrifuge Tube with flat cap, Amber	Bulk	25 Tubes	25 Tubes X 20 packs
ACT-50-K-B-S-A	50 ml Self Standing Centrifuge Tube with flat Cap, Sterile, Amber	Bulk	25 Tubes	25 Tubes X 20 packs
ACT-50-K-S-A-W	50 ml Self Standing Centrifuge Tube with flat Cap, Sterile, Amber, Individually Wrapped	Bulk	25 Tubes	25 Tubes X 20 packs
ACT-50-K-B	50 ML Self Standing Centrifuge Tube With Flat Cap	Bulk	25 Tubes	25 Tubes X 20 packs
ACT-50-K-B-S	50 ML Self Standing Centrifuge Tube With Flat Cap, Sterile	Bulk	25 Tubes	25 Tubes X 20 packs
ACT-50-K-S-W	50 ML Self Standing Centrifuge Tube With Flat Cap, Sterile, Individually Wrapped	Individually Wrapped	25 Tubes	25 Tubes x 20 packs
ACT-50-K-R-A	50 ml Self Standing Centrifuge Tube with flat cap, Amber	Rack	25 Tubes	25 Tubes X 20 packs
ACT-50-K-R-S-A	50 ml Self Standing Centrifuge Tube with flat cap, Sterile Amber	Rack	25 Tubes	25 Tubes X 20 packs

## Accumax Metal Free Tubes

Cat. No.	Description	Packing Type	Quantity / Pack	Quantity / Case
ACT-15-B-T	15 ML Centrifuge Tube With Heavy Metal Free Flat Cap	Bulk	25 Tubes	25 Tubes X 20 packs
ACT-15-B-S-T	15 ML Centrifuge Tube With Heavy Metal Free Flat Cap, Sterile	Rack	25 Tubes	25 Tubes X 20 packs
ACT-50-B-T	50 ML Centrifuge Tube With Heavy Metal Free Flat Cap	Bulk	25 Tubes	25 Tubes X 20 packs
ACT-50-B-S-T	50 ML Centrifuge Tube With Heavy Metal Free Flat Cap, Sterile	Rack	25 Tubes	25 Tubes X 20 packs
ACT-50-K-B-T	50ml Self Standing Centrifuge Tube With Heavy Metal Free Flat Cap	Bulk	25 Tubes	25 Tubes X 20 packs
ACT-50-K-B-S-T	50ml Self Standing Centrifuge Tube With Heavy Metal Free Flat Cap , Sterile	Bulk	25 Tubes	25 Tubes X 20 packs
ACT-50-K-S-W-T	50ml Self Standing Centrifuge Tube With Heavy Metal Free Flat Cap, Sterile, Individually Wrapped	Individually Wrapped	25 Tubes	25 Tubes X 20 packs

# ACCUMAX

**Accumax Lab Devices Pvt. Ltd.**  
www.accumaximum.com

ISO 9001:2015  
ISO 13485:2016

\*We also execute OEM projects. Specifications can be changed without notice for quality improvement.