


ACCUMAX




MAJESTIC COMPATIBILITY
DEEPWELL PLATES

Accumax Deepwell Plates

Accumax Deepwell Plates

In high-precision environment sample storage has to be as efficient as the instruments of measurement and many manufacturers today fail to pay attention to these miniscule requirements. As a result, chances of sample loss and measurement errors increase.


With Accumax Deepwell Plates you can forget about inefficient sample measurement or recovery as our cutting edge product is designed after thorough R&D. Our unique features like impeccable temperature tolerance or alphanumeric indexing minimizes any chance of error in measurement or in sample recovery.




**2.2 ML, ZYBIO®
TYPE U-BOTTOM**




**1 ML ROUND
BOTTOM**



2 ML V-BOTTOM



**2.2 ML, KINGFISHER®
TYPE V-BOTTOM**



**2 ML ROUND
BOTTOM**



1 ML U-Bottom

*Kingfisher® Kingfisher is registered trademark owned by Thermo Scientific, USA & other countries

*Zybio® Zybio is registered owned by Zybio Inc., China & other countries

Exceptional Excellence



Manufactured In
A Fully Automated
Class 100K Clean Room




Ensures Chemical Resistance
And Sample Integrity With
Long-term Storage Needs.



Free Of Detectable
DNase, RNase, And
Endotoxin Inhibitors.




Comes With An Easy
And Reliable Stacking.



Raised rim design to
avoid contamination


Fully Autoclavable



Provides Temperature
Resistance Up To -80°C




Alphanumeric Indexing
For Quick And Easy
Sample Identification.




Made Considering All
Important ANSI Standards.




Cross-compatible silicon mat available in square and round design



Silicon mat all fit for
square top cross cut



Silicon mat round
1.0 ml pierceable



Silicon mat round
2.0 ml pierceable

Deepwell Compatibility Matrix

Deepwell Compatibility Matrix							
Accumax Cat Code	Deepwell Design Types	Thermo kingfisher APEX	Thermo kingfisher FLEX presto 32 96	Thermo kingfisher presto 24/ 96	Allsheng auto PURE-16	Allsheng auto PURE-32	Allsheng auto PURE-96
ADW-2V	2.0 ml, V-bottom-Square top deep well plate						
ADW-1R	1.0 ml, Round bottom-round top deep well plate						
ADW-2.2-U	2.2 ml , U-bottom zybio type - Square top				☑	☑	
ADW-TC-08	8 well tip comb - for 96 U bottom zybio type 2.2 ml deep well plate						
ADW-2.2-K	2.2 ml, Kingfisher type bottom - Square top deep well plate	☑	☑	☑			☑
ADW-TC-96	96 well Kingfisher type tip comb - for 2.2 ml plate	☑	☑	☑			☑
ADW-2R	2.0 ml, U-bottom-round top deep well plate						
ADW-1U	1.0 ml, U-bottom Square top deep well plate						

Deepwell Compatibility Matrix							
Accumax Cat Code	Deepwell Design Types	Zybio EXM3000 32	Zybio EXM6000 96	Himedia InstaNX Mag 32	Himedia InstaNX Mag 96	Q line Biotech POCT services Q Xtra 32	Q line Biotech POCT services Q Xtra 96
ADW-2V	2.0 ml, V-bottom-Square top deep well plate						
ADW-1R	1.0 ml, Round bottom-round top deep well plate						
ADW-2.2-U	2.2 ml , U-bottom zybio type - Square top			☑			
ADW-TC-08	8 well tip comb - for 96 U bottom zybio type 2.2 ml deep well plate	☑				☑	
ADW-2.2-K	2.2 ml, Kingfisher type bottom - Square top deep well plate		☑		☑		☑
ADW-TC-96	96 well Kingfisher type tip comb - for 2.2 ml plate		☑		☑		☑
ADW-2R	2.0 ml, U-bottom-round top deep well plate						
ADW-1U	1.0 ml, U-bottom Square top deep well plate						

Deepwell Compatibility Matrix							
Accumax Cat Code	Deepwell Design Types	Labgene ideal 32	Labgene ideal 96	Trangene biotech TS 32	Trangene biotech TS 96	Genetix Purifier 32	Genetix Purifier HT 96
ADW-2V	2.0 ml, V-bottom-Square top deep well plate						
ADW-1R	1.0 ml, Round bottom-round top deep well plate						
ADW-2.2-U	2.2 ml , U-bottom zybio type - Square top	☑	☑			☑	
ADW-TC-08	8 well tip comb - for 96 U bottom zybio type 2.2 ml deep well plate			☑			
ADW-2.2-K	2.2 ml, Kingfisher type bottom - Square top deep well plate				☑		☑
ADW-TC-96	96 well Kingfisher type tip comb - for 2.2 ml plate				☑		☑
ADW-2R	2.0 ml, U-bottom-round top deep well plate						
ADW-1U	1.0 ml, U-bottom Square top deep well plate						

Deepwell Compatibility Matrix

Deepwell Compatibility Matrix							
Accumax Cat Code	Deepwell Design Types	id solution IDEAL-32	id solution IDEAL-96	Bioevopeak NAE-9032	Bioevopeak NAE-0196	BIOER Gene Pure Plus 32	BIOER Gene Pure Pro 96
ADW-2V	2.0 ml, V-bottom-Square top deep well plate						
ADW-1R	1.0 ml, Round bottom-round top deep well plate						
ADW-2.2-U	2.2 ml , U-bottom zybio type - Square top	☑	☑	☑		☑	
ADW-TC-08	8 well tip comb - for 96 U bottom zybio type 2.2 ml deep well plate					☑	
ADW-2.2-K	2.2 ml, Kingfisher type bottom - Square top deep well plate				☑		☑
ADW-TC-96	96 well Kingfisher type tip comb - for 2.2 ml plate				☑		
ADW-2R	2.0 ml, U-bottom-round top deep well plate						
ADW-1U	1.0 ml, U-bottom Square top deep well plate						

Deepwell Compatibility Matrix							
Accumax Cat Code	Deepwell Design Types	MGI SP-NE32	MGISP-NE384	Dianbio EB-1000 32	Biobase BNP 32	Biobase HS 96	Genosolution Nexttractor® NX-48S
ADW-2V	2.0 ml, V-bottom-Square top deep well plate						
ADW-1R	1.0 ml, Round bottom-round top deep well plate						
ADW-2.2-U	2.2 ml , U-bottom zybio type - Square top	☑					
ADW-TC-08	8 well tip comb - for 96 U bottom zybio type 2.2 ml deep well plate	☑			☑		☑
ADW-2.2-K	2.2 ml, Kingfisher type bottom - Square top deep well plate	☑	☑			☑	
ADW-TC-96	96 well Kingfisher type tip comb - for 2.2 ml plate		☑			☑	
ADW-2R	2.0 ml, U-bottom-round top deep well plate						
ADW-1U	1.0 ml, U-bottom Square top deep well plate						

Deepwell Compatibility Matrix							
Accumax Cat Code	Deepwell Design Types	Zhuai HEMA E96-II 96	Hyribio HBNP-4801A 48	Tiangen YOSE-S32 32	Tiangen OSE-S96-01 96	Bioteke AU1001 32	Bioteke AU1001-96
ADW-2V	2.0 ml, V-bottom-Square top deep well plate						
ADW-1R	1.0 ml, Round bottom-round top deep well plate						
ADW-2.2-U	2.2 ml , U-bottom zybio type - Square top		☑	☑		☑	☑
ADW-TC-08	8 well tip comb - for 96 U bottom zybio type 2.2 ml deep well plate	☑					
ADW-2.2-K	2.2 ml, Kingfisher type bottom - Square top deep well plate				☑		
ADW-TC-96	96 well Kingfisher type tip comb - for 2.2 ml plate						
ADW-2R	2.0 ml, U-bottom-round top deep well plate						
ADW-1U	1.0 ml, U-bottom Square top deep well plate						

Deepwell Compatibility Matrix

Deepwell Compatibility Matrix							
Accumax Cat Code	Deepwell Design Types	Backed SSNP-2000A 32	XATNKJ NP-968 32	Zeean Lad-Aid 808 8	Abclonal ANDIS 350	Jiyuan biotech BFEX-32	Geneunion BFEX-96
ADW-2V	2.0 ml, V-bottom-Square top deep well plate						
ADW-1R	1.0 ml, Round bottom-round top deep well plate						
ADW-2.2-U	2.2 ml , U-bottom zybio type - Square top	☑		☑	☑	☑	
ADW-TC-08	8 well tip comb - for 96 U bottom zybio type 2.2 ml deep well plate		☑		☑		
ADW-2.2-K	2.2 ml, Kingfisher type bottom - Square top deep well plate						☑
ADW-TC-96	96 well Kingfisher type tip comb - for 2.2 ml plate						☑
ADW-2R	2.0 ml, U-bottom-round top deep well plate						
ADW-1U	1.0 ml, U-bottom Square top deep well plate						

Deepwell Compatibility Matrix							
Accumax Cat Code	Deepwell Design Types	Denogen DOF-9696PRO	Bioneer Exiprep 16 plus	Bioneer Exiprep 96 lite	Bioenriching ETP96-X2 beads ETP-300	TIAGEN 192 SAMPLE S-96	Lepu medical Nexor 32
ADW-2V	2.0 ml, V-bottom-Square top deep well plate						
ADW-1R	1.0 ml, Round bottom-round top deep well plate						
ADW-2.2-U	2.2 ml , U-bottom zybio type - Square top		☑	☑	☑		☑
ADW-TC-08	8 well tip comb - for 96 U bottom zybio type 2.2 ml deep well plate						
ADW-2.2-K	2.2 ml, Kingfisher type bottom - Square top deep well plate	☑				☑	
ADW-TC-96	96 well Kingfisher type tip comb - for 2.2 ml plate						
ADW-2R	2.0 ml, U-bottom-round top deep well plate						
ADW-1U	1.0 ml, U-bottom Square top deep well plate						

Deepwell Compatibility Matrix							
Accumax Cat Code	Deepwell Design Types	Lepu medical Nexor 96	Esco Swift Ex32	Lovley perfume store ZK01 32	Lovley perfume store Zk96	MN Filtration chromabond multi	Qiagen Filtration qiavac 96
ADW-2V	2.0 ml, V-bottom-Square top deep well plate					☑	☑
ADW-1R	1.0 ml, Round bottom-round top deep well plate						
ADW-2.2-U	2.2 ml , U-bottom zybio type - Square top		☑				
ADW-TC-08	8 well tip comb - for 96 U bottom zybio type 2.2 ml deep well plate		☑	☑			
ADW-2.2-K	2.2 ml, Kingfisher type bottom - Square top deep well plate	☑				☑	
ADW-TC-96	96 well Kingfisher type tip comb - for 2.2 ml plate	☑				☑	
ADW-2R	2.0 ml, U-bottom-round top deep well plate						
ADW-1U	1.0 ml, U-bottom Square top deep well plate						

Deepwell Compatibility Matrix

Deepwell Compatibility Matrix							
Accumax Cat Code	Deepwell Design Types	HPLC like Agilent Infinity sample prep	LC MS like Shimadzu sample prep	GC MS like perkin clarus sq 8c sample prep	Ultravap mistral systems	Waters plate manifold oasis & other systems	Thermo Cryomed
ADW-2V	2.0 ml, V-bottom-Square top deep well plate	✓	✓	✓	✓	✓	✓
ADW-1R	1.0 ml, Round bottom-round top deep well plate	✓	✓	✓	✓	✓	✓
ADW-2.2-U	2.2 ml , U-bottom zybio type - Square top						
ADW-TC-08	8 well tip comb - for 96 U bottom zybio type 2.2 ml deep well plate						
ADW-2.2-K	2.2 ml, Kingfisher type bottom - Square top deep well plate						
ADW-TC-96	96 well Kingfisher type tip comb - for 2.2 ml plate						
ADW-2R	2.0 ml, U-bottom-round top deep well plate	✓	✓	✓	✓	✓	✓
ADW-1U	1.0 ml, U-bottom Square top deep well plate	✓	✓	✓	✓	✓	✓

Deepwell Compatibility Matrix							
Accumax Cat Code	Deepwell Design Types	Thermo Cryoplus	Epigenetik magnetic EPI Mag HT	Roche Magnapure 24	Roche Magnapure 96	Eppendorf Ep Motion	Caliper perkinelmer Janus
ADW-2V	2.0 ml, V-bottom-Square top deep well plate	✓		✓	✓	✓	✓
ADW-1R	1.0 ml, Round bottom-round top deep well plate	✓	✓	✓	✓	✓	✓
ADW-2.2-U	2.2 ml , U-bottom zybio type - Square top						
ADW-TC-08	8 well tip comb - for 96 U bottom zybio type 2.2 ml deep well plate						
ADW-2.2-K	2.2 ml, Kingfisher type bottom - Square top deep well plate						
ADW-TC-96	96 well Kingfisher type tip comb - for 2.2 ml plate						
ADW-2R	2.0 ml, U-bottom-round top deep well plate	✓	✓	✓	✓	✓	✓
ADW-1U	1.0 ml, U-bottom Square top deep well plate	✓	✓	✓	✓	✓	✓

Deepwell Compatibility Matrix							
Accumax Cat Code	Deepwell Design Types	Caliper perkinelmer Sciclone G3	Caliper perkinelmer Zephyr G3	Perkinelmer Chemagic 360	Tecan Fluent	Tecan Freedom Evo	Beckman Biomek i5
ADW-2V	2.0 ml, V-bottom-Square top deep well plate	✓	✓	✓	✓	✓	✓
ADW-1R	1.0 ml, Round bottom-round top deep well plate	✓	✓	✓	✓	✓	✓
ADW-2.2-U	2.2 ml , U-bottom zybio type - Square top						
ADW-TC-08	8 well tip comb - for 96 U bottom zybio type 2.2 ml deep well plate						
ADW-2.2-K	2.2 ml, Kingfisher type bottom - Square top deep well plate						
ADW-TC-96	96 well Kingfisher type tip comb - for 2.2 ml plate						
ADW-2R	2.0 ml, U-bottom-round top deep well plate	✓	✓	✓	✓	✓	✓
ADW-1U	1.0 ml, U-bottom Square top deep well plate	✓	✓	✓	✓	✓	✓

Deepwell Compatibility Matrix

Deepwell Compatibility Matrix								
Accumax Cat Code	Deepwell Design Types	Beckman Biomek i7	Beckman Biomek 4000	Biomek NGenius	Opentron	Spexs Agriprep homogenizers genogrinders	Dr. oligo Biolytic	Hamilton based Liquid handlers
ADW-2V	2.0 ml, V-bottom-Square top deep well plate	✓	✓	✓	✓	✓	✓	✓
ADW-1R	1.0 ml, Round bottom-round top deep well plate	✓	✓	✓	✓	✓	✓	✓
ADW-2.2-U	2.2 ml , U-bottom zybio type - Square top							
ADW-TC-08	8 well tip comb - for 96 U bottom zybio type 2.2 ml deep well plate							
ADW-2.2-K	2.2 ml, Kingfisher type bottom - Square top deep well plate							
ADW-TC-96	96 well Kingfisher type tip comb - for 2.2 ml plate							
ADW-2R	2.0 ml, U-bottom-round top deep well plate	✓	✓	✓	✓	✓	✓	✓
ADW-1U	1.0 ml, U-bottom Square top deep well plate	✓	✓	✓	✓	✓	✓	✓

Ordering Information

Accumax Deepwell Plates			
Cat. No.	Description	Unit/Pack	Unit/Case
ADW-2V	2.0 ml, V-bottom deep well plate, Bulk	12/Pack	48/Case
ADW-1R	1.0 ml, Round-bottom deep well plate, Bulk	24/Pack	96/Case
ADW-2V-S	2.0 ml, Sterile Bulk	12/Pack	48/Case
ADW-1R-S	1.0 ml, Sterile Bulk	24/Pack	96/Case
ADW-2V-S-W	2.0 ml, Sterile individual wrap	12/Pack	48/Case
ADW1R-S-W	1.0 ml, Sterile individual wrap	24/Pack	96/Case
ADW-2R	2.0 ml, U-bottom deep well plate	12/Pack	48/Case
ADW-2R-S	2.0 ml , U-bottom deep well plate, Sterile, Bulk	12/Pack	48/Case
ADW-2.2 K	2.2 ml, Kingfisher® type V bottom - Square top deep well plate, bulk	5/Pack	10/case, 50 No.
ADW-2.2 U	2.2 ml, zybio® type U-bottom - Square top deep well plate, bulk	5/Pack	10/case, 50 No.
ADW-TC-96	Kingfisher® type tip comb - for 96 well 2.2 ml plate	2/Pack	25/case, 50 No.
ADW-TC-08	8 well tip comb - for 96 U bottom zybio® type 2.2 ml deep well plate	2/Pack	100/case, 200 No.
ADW-1.0 U	1.0 ml, U bottom Square top deep well plate, bulk	24/Pack	4 Packs/case
ADW-1.0 U-S	1.0 ml, U bottom Square top deep well plate, Sterile, bulk	24/Pack	4 Packs/case
ADW-1.0 U-S-W	1.0 ml, U bottom Square top deep well plate, sterile, individually wrapped	24/Pack	4 Packs/case

*Kingfisher® Kingfisher is registered trademark owned by Thermo Scientific, USA & other countries

*Zybio® Zybio is registered owned by Zybio Inc., China & other countries

Ordering Information

Accumax Deepwell Plates			
Cat. No.	Description	Unit/Pack	Unit/Case
ADW-SM-CC	Silicon mat all fit for square top cross cut	10/Pack	50/Case
ADW-SM-1-RP	Silicon mat round 1.0 ml pierceable	10/Pack	50/Case
ADW-SM-2-RP	Silicon mat round 2.0 ml pierceable	10/Pack	50/Case

Accumax Deepwell Dilution Strips		
Cat. No.	Description	Unit/Pack/Case
ADTD-2	Deepwell Strips 2ml	250

ACCUMAX

Accumax Lab Devices Pvt. Ltd.
www.accumaxlab.com

ISO 9001:2015
ISO 13485:2016

*We also execute OEM projects. Specifications can be changed without notice for quality improvement.