

FILTRATION RANGE

Microfilter Plates

Accumax




Sample volume of 1 ml



Certified free from RNase/DNase, human DNA, pyrogens & PCR inhibitors with sterility assurance level (SAL) of 10^{-6}



Manufactured in ISO 9001 & ISO 13485 certified 100K clean room facility



Crafted from choice-grade virgin polystyrene, ensuring each microfilter plate has outstanding transparency



Multiple filter options like Glass Fibre filter, Polyethylene filter, Cellulose Acetate filter, Regenerated Cellulose, Nylon filter, Polyvinylidene Fluoride filter




Available in brilliant options for prefiltration, coarse filtration, DNA isolation, RNA isolation, protein purification, PCR clean up and chromatography sample preparations



Compatible with centrifugation (1500g) and vacuum protocols (10 inHG)



Fitted with long drip directors that allow for secure use while concomitantly reducing cross-contamination



For more details visit us at www.accumaximum.com