

ACCUMAX




LOW FORCE AND MAXIMUM
EFFICIENCY


ACCUMAX LF TIPS

ACCUMAX LF TIPS

Accumax LF pipette tips known for Low Force with maximum efficiency are designed to fit Rainin® LTS® style pipette, providing accuracy, precision and ergonomics you have come to expect when using your RAININ® style pipette. It provides a perfect seal between the pipette and tips with minimal effort.



Cylindrical shaped for the perfect fit with minimal effort



PRODUCT FEATURES



Designed for optimal performance with Rainin® LTS® style pipette




Certified Free of RNase, DNase, PCR Inhibitors, Endotoxins, and Pyrogens




100k clean room facility with automated manufacturing

Cylindrical shaped and positive-stop design ensures a perfect fit,



Color-coded inserts for easy matching of the tip



US FDA approved medical grade virgin polypropylene plastic without any trace of metals, highly transparent and clear



PRODUCT FEATURES



All tips have graduation marks which act as visual checks to minimize any risk of errors during sample aspiration




Rack, Bulk, Reload and ECO rack






Swiss-made diamond finish mold for flawless "super smooth" internal surface



Specially developed PE filters stop aerosols contamination from entering into pipette and thus lessens the risk of cross-contamination








Ordering Information

Cat#	Description	Tray color	Type of packing	Quantity/case	Quantity / Pack(1 Rack =96 tips)	Qty/Case
ALF-20-R	20µL Tips for Low Force (LF)Pipettes	Red	Rack	96 Tips x 10 Racks x 10 Packs	10 Racks	10 Packs
ALF-20-S-R	20µL Sterile Tips for Low Force (LF)Pipettes	Red	Rack	96 Tips x 10 Racks x 10 Packs	10 Racks	10 Packs
ALF-20-S-L-R	20µL Sterile Low-Retention Tips for Low Force (LF)Pipettes	Red	Rack	96 Tips x 10 Racks x 10 Packs	10 Racks	10 Packs
ALF-20-S-F-R	20µL Sterile Filter Tips for Low Force (LF)Pipettes	Red	Rack	96 Tips x 10 Racks x 10 Packs	10 Racks	10 Packs
ALF-20-S-L-F-R	20µL Sterile Low-Retention Filter Tips for Low Force (LF)Pipettes	Red	Rack	96 Tips x 10 Racks x 10 Packs	10 Racks	10 Packs
ALF-20-K	20µL Tips for Low Force (LF)Pipettes	Red	Reloading stack	96 Tips x 10 Racks x 10 Packs	10 Racks	10 Packs
ALF-20-S-K	20µL Sterile Tips for Low Force (LF)Pipettes	Red	Reloading stack	96 Tips x 10 Racks x 10 Packs	10 Racks	10 Packs
ALF-20-S-L-K	20µL Sterile Low-Retention Tips for Low Force (LF)Pipettes	Red	Reloading stack	96 Tips x 10 Racks x 10 Packs	10 Racks	10 Packs
GLF-20-R	20µL Tips for Low Force (LF)Pipettes 	Red	Eco-Rack	96 Tips x 10 Racks x 10 Packs	10 Racks	10 Packs
GLF-20-S-R	20µL Sterile Tips for Low Force (LF)Pipettes 	Red	Eco-Rack	96 Tips x 10 Racks x 10 Packs	10 Racks	10 Packs
GLF-20-S-F-R	20µL Sterile Filter Tips for Low Force (LF)Pipettes 	Red	Eco-Rack	96 Tips x 10 Racks x 10 Packs	10 Racks	10 Packs
GLF-20-S-L-R	20µL Sterile Low-Retention Tips for Low Force (LF)Pipettes 	Red	Eco-Rack	96 Tips x 10 Racks x 10 Packs	10 Racks	10 Packs
GLF-20-S-L-F-R	20µL Sterile Low-Retention Filter Tips for Low Force (LF)Pipettes 	Red	Eco-Rack	96 Tips x 10 Racks x 10 Packs	10 Racks	10 Packs
ALF-20-B	20µL Bulk tips for Low Force (LF)Pipette	NA	Bulk	1000 tips x 40 packs	1000 tips	40 Packs
ALF-20-S-B	20µL Sterile Bulk tips for Low Force (LF)Pipette	NA	Bulk	1000 tips x 40 packs	1000 tips	40 Packs
ALF-20-F-B	20µL Bulk Filter tips for Low Force (LF)Pipette	NA	Bulk	1000 tips x 40 packs	1000 tips	40 Packs
ALF-20-S-F-B	20µL Bulk Sterile Filter tips for Low Force (LF)Pipette	NA	Bulk	1000 tips x 40 packs	1000 tips	40 Packs
ALF-20-S-L-B	20µL Sterile Low-Retention Bulk tips for Low Force (LF)Pipette	NA	Bulk	1000 tips x 40 packs	1000 tips	40 Packs

Compatible to Rainin® LTS® pipettes

LTS® is a registered trademark of Rainin Instrument






Rainin® is a registered trademark of Mettler Toledo Rainin LLC

Cat#	Description	Tray color	Type of packing	Quantity/case	Quantity / Pack(1 Rack =96 tips)	Qty/Case
ALF-20-S-L-F-B	20µL Sterile Low-Retention Filter Bulk tips for Low Force (LF)Pipette	NA	Bulk	1000 tips x 40 Packs	1000 tips	40 Packs
ALF-200-R	200µL Tips for Low Force (LF)Pipettes	Green	Rack	96 Tips x 10 Racks x 10 Packs	10 Racks	10 Packs
ALF-200-S-R	200µL Sterile Tips for Low Force (LF)Pipettes	Green	Rack	96 Tips x 10 Racks x 10 Packs	10 Racks	10 Packs
ALF-200-S-L-R	200µL Sterile Low-Retention Tips for Low Force (LF)Pipettes	Green	Rack	96 Tips x 10 Racks x 10 Packs	10 Racks	10 Packs
ALF-200-S-F-R	200µL Sterile Filter Tips for Low Force (LF)Pipettes	Green	Rack	96 Tips x 10 Racks x 10 Packs	10 Racks	10 Packs
ALF-200-S-L-F-R	200µL Sterile Low-Retention Filter Tips for Low Force (LF)Pipettes	Green	Rack	96 Tips x 10 Racks x 10 Packs	10 Racks	10 Packs
ALF-200-K	200µL Tips for Low Force (LF)Pipettes	Green	Reloading stack	96 Tips x 10 Racks x 10 Packs	10 Racks	10 Packs
ALF-200-S-K	200µL Sterile Tips for Low Force (LF)Pipettes	Green	Reloading stack	96 Tips x 10 Racks x 10 Packs	10 Racks	10 Packs
ALF-200-S-L-K	200µL Sterile Low-Retention Tips for Low Force (LF)Pipettes	Green	Reloading stack	96 Tips x 10 Racks x 10 Packs	10 Racks	10 Packs
GLF-200-R	200µL Tips for Low Force (LF)Pipettes 	Green	Eco-Rack	96 Tips x 10 Racks x 10 Packs	10 Racks	10 Packs
GLF-200-S-R	200µL Sterile Tips for Low Force (LF)Pipettes 	Green	Eco-Rack	96 Tips x 10 Racks x 10 Packs	10 Racks	10 Packs
GLF-200-S-F-R	200µL Sterile Filter Tips for Low Force (LF)Pipettes 	Green	Eco-Rack	96 Tips x 10 Racks x 10 Packs	10 Racks	10 Packs
GLF-200-S-L-R	200µL Sterile Low retention Tips for Low Force (LF)Pipettes 	Green	Eco-Rack	96 Tips x 10 Racks x 10 Packs	10 Racks	10 Packs
GLF-200-S-L-F-R	200µL Sterile Low retention Filter Tips for Low Force (LF)Pipettes 	Green	Eco-Rack	96 Tips x 10 Racks x 10 Packs	10 Racks	10 Packs
ALF-200-B	200µL Bulk Tips for Low Force (LF)Pipettes	NA	Bulk	1000 tips x 20 packs	1000 tips	20 Packs
ALF-200-S-B	200µL Sterile Bulk Tips for Low Force (LF)Pipettes	NA	Bulk	1000 tips x 20 packs	1000 tips	20 Packs
ALF-200-F-B	200µL Filter Bulk tips for Low Force (LF)Pipette	NA	Bulk	1000 tips x 20 packs	1000 tips	20 Packs
ALF-200-S-F-B	200µL Sterile Filter Bulk tips for Low Force (LF)Pipette	NA	Bulk	1000 tips x 20 packs	1000 tips	20 Packs

Compatible to Rainin® LTS® pipettes

LTS® is a registered trademark of Rainin Instrument

Rainin® is a registered trademark of Mettler Toledo Rainin LLC

Cat#	Description	Tray color	Type of packing	Quantity/case	Quantity / Pack(1 Rack =96 tips)	Qty/Case
ALF-200-S-L-B	200µL Sterile Low retention Bulk tips for Low Force (LF)Pipette	NA	Bulk	1000 tips x 20 Packs	1000 tips	20 Packs
ALF-200-S-L-F-B	200µL Sterile Low retention Filter Bulk tips for Low Force (LF)Pipette	NA	Bulk	1000 tips x 40 Packs	1000 tips	20 Packs
ALF-1000-R	1000µL Tips for Low Force (LF)Pipettes	Blue	Rack	96 Tips x 10 Racks x 10 Packs	6 Racks	10 Packs
ALF-1000-S-R	1000µL Sterile Tips for Low Force (LF)Pipettes	Blue	Rack	96 Tips x 10 Racks x 10 Packs	6 Racks	10 Packs
ALF-1000-S-L-R	1000µL Sterile Low-Retention Tips for Low Force (LF)Pipettes	Blue	Rack	96 Tips x 10 Racks x 10 Packs	6 Racks	10 Packs
ALF-1000-S-F-R	1000µL Sterile Filter Tips for Low Force (LF)Pipettes	Blue	Rack	96 Tips x 10 Racks x 10 Packs	6 Racks	10 Packs
ALF-1000-S-L-F-R	1000µL Sterile Low-Retention Filter Tips for Low Force (LF)Pipettes	Blue	Rack	96 Tips x 10 Racks x 10 Packs	6 Racks	10 Packs
ALF-1000-K	1000µL Tips for Low Force (LF)Pipettes	Blue	Reloading stack	96 Tips x 10 Racks x 10 Packs	6 Racks	10 Packs
ALF-1000-S-K	1000µL Sterile Tips for Low Force (LF)Pipettes	Blue	Reloading stack	96 Tips x 10 Racks x 10 Packs	6 Racks	10 Packs
ALF-1000-S-L-K	1000µL Sterile Low-Retention Tips for Low Force (LF)Pipettes	Blue	Reloading stack	96 Tips x 10 Racks x 10 Packs	6 Racks	10 Packs
GLF-1000-R	1000µL Tips for Low Force (LF)Pipettes 	Blue	Eco-Rack	96 Tips x 10 Racks x 10 Packs	6 Racks	10 Packs
GLF-1000-S-R	1000µL Sterile Tips for Low Force (LF)Pipettes 	Blue	Eco-Rack	96 Tips x 10 Racks x 10 Packs	6 Racks	10 Packs
GLF-1000-S-F-R	1000µL Sterile Filter Tips for Low Force (LF)Pipettes 	Blue	Eco-Rack	96 Tips x 10 Racks x 10 Packs	6 Racks	10 Packs
GLF-1000-S-L-R	1000µL Sterile Low-Retention Tips or Low Force (LF)Pipettes 	Blue	Eco-Rack	96 Tips x 10 Racks x 10 Packs	6 Racks	10 Packs
GLF-1000-S-L-F-R	1000µL Sterile Low-Retention Filter Tips for Low Force (LF)Pipettes 	Blue	Eco-Rack	96 Tips x 10 Racks x 10 Packs	6 Racks	10 Packs
ALF-1000-B	1000µL Bulk Tips for Low Force (LF)Pipettes	NA	Bulk	1000 tips x 40 packs	1000 tips	10 Packs
ALF-1000-S-B	1000µL Sterile Bulk Tips for Low Force (LF)Pipettes	NA	Bulk	1000 tips x 40 packs	1000 tips	10 Packs
ALF-1000-F-B	1000µL Filter Bulk Tips for Low Force (LF)Pipettes	NA	Bulk	1000 tips x 40 packs	1000 tips	10 Packs

Compatible to Rainin® LTS® pipettes

LTS® is a registered trademark of Rainin Instrument

Rainin® is a registered trademark of Mettler Toledo Rainin LLC

Cat#	Description	Tray color	Type of packing	Quantity/case	Quantity / Pack(1 Rack =96 tips)	Qty/Case
ALF-1000-S-F-B	1000µL Sterile Filter Bulk Tips for Low Force (LF)Pipettes	NA	Bulk	1000 tips x 40 packs	1000 tips	10 Packs
ALF-1000-S-L-B	1000µL Sterile Low retention Bulk Tips for Low Force (LF)Pipettes	NA	Bulk	1000 tips x 10 Packs	1000 tips	10 Packs
ALF-1000-S-L-F-B	1000µL Sterile Low retention Filter Bulk Tips for Low Force (LF)Pipettes	NA	Bulk	1000 tips x 10 Packs	1000 tips	10 Packs
ALF-20-E-R	20µL Empty rack tip - Red color tray for LTS pipette with out tips	NA	Bulk	10 Racks x 10 Packs	10 Racks	10 Packs
ALF-200-E-R	200µL Empty rack tip - Green color tray for LTS pipette with out tips	NA	Bulk	10 Racks x 10 Packs	10 Racks	10 Packs
ALF-1000-E-R	1000µL Empty rack tip - Blue color tray for LTS pipette with out tips	NA	Bulk	10 Racks x 10 Packs	10 Racks	10 Packs

Compatible to Rainin® LTS® pipettes

LTS® is a registered trademark of Rainin Instrument

Rainin® is a registered trademark of Mettler Toledo Rainin LLC

ACCUMAX

Accumax Lab Devices Pvt. Ltd.
www.accumaximum.com

ISO 9001:2015
ISO 13485:2016

*We also execute OEM projects. Specifications can be changed without notice for quality improvement.