

VERSATILE, RELIABLE AND CONSISTENT

ACCUMAX

TECAN® TYPE TIPS

Volume: 10,50,200 and 1000 µl




Constructed of premium-grade virgin polypropylene. Available in carbon-embedded conductive tips and clear tips




USP class VI resins



Validated on the actual instrument for which the tip is created



Polycarbonate (PC) tray with color coding to indicate different volumes



Available in ANSI/SLAS Rack, Blister packs & Eco flip box




Manufactured in ISO- certified 9001 and 13485 facility







All tips are certified Endotoxin, Heavy Metals, RNase, DNase and DNA free



Nonsterile, sterile and filtered sterile formats



Polyethylene filter to prevent aerosol contamination



TECAN® is registered trademark of Tecan Trading AG

For more details visit us at www.accumaxlab.com

AN ISO 9001 : 2015 & ISO 13485 : 2016 CERTIFIED MFG. COMPANY

*We also execute OEM projects. Specificness can be changed without notice for quality improvement.