

SECOND TO NONE DEEPWELL PLATES

Accumax Deepwell Plates

ACCUMAX




Versatile 96 well format for small volume biological sample operation and storage



Alphanumeric pattern for simple identification of sample



Overall height of 33mm and 44mm respectively for 1.2 & 2.2ml DWP



Fully autoclavable and temperature resistance from -80°C to 121°C



Fits to hamilton liquid handling automation machine




Designed as per SBS and ANSI standards for broad instrument compatibility




U and V bottom for optimal sample recovery



1.1 ml round bottom deepwell plate



2.2 ml V-bottom deepwell plate



Hamilton star using deepwell plate