

ACCUMAX

**THE COMPACT
POWERHOUSE**

Accumax Jr. Pipette

Accumax Jr. Pipette

Research indicates that laboratories are in need of a compact, quick and effortless pipetting device that can help them perform small tasks with utter ease. Moreover even small laboratories and college labs cannot afford to buy specialised equipment for training & experimenting purposes.

The Accumax Jr. Pipette is a device born out of that need. Accumax takes a holistic approach to product development by engaging in extensive research about user data which informs the design of the pipette. The result are pipettes that deliver exceptional performance in a compact package with great ergonomics for a great user experience.



Product Features




Compact Working Size Of 130 Mm Ensures Fatigue Free Operation



Two Step Plunger Operation Allows For "Reverse Pipetting Technique



Calibrated As Per En Iso 8655-6 Standard To Ensure Relentless Precision & Accuracy



Fully Autoclavable Body Makes It Easy To Sterilize

Specially Designed Tip Cone For Both - 10 μ l & 200 μ l Tips To Ensure Convenient & Effortless Usage.



Stainless Steel Piston For Reliable Performance And Years Of Durability.



Endless Customization In Terms Of Any Volume Between 2 μ l To 100 μ l And Color Combination For Supply Under Private Label



Available In Different Volume Range From 5 μ l To 200 μ l



Ordering Information

Accumax Jr. Pipette				
Cat. No.	Volume Range (µL)	Test Volume	Inaccuracy (±) %	Imprecision (±) %
AJ-5	5	5	1.3	1.2
AJ-10	10	10	0.8	0.8
AJ-20	20	20	0.6	0.5
AJ-25	25	25	0.5	0.3
AJ-50	50	50	0.5	0.3
AJ-100	100	100	0.5	0.3
AJ-200	200	200	0.4	0.2

ACCUMAX

Accumax Lab Devices Pvt. Ltd.
www.accumaximum.com

ISO 9001:2015
ISO 13485:2016

*We also execute OEM projects. Specifications can be changed without notice for quality improvement.