

ROBUST RESISTANT & RELIABLE

Reagent Bottles (PP* /HDPE/LDPE)

Narrow Mouth: 4/8/15/30/60/125/250/500ml

Wide Mouth: 30/60/125/250/500ml

ACCUMAX



Available in wide and narrow Mouth, in different volumes

Amber colour for excellent light shielding for storage of light sensitive substance



Excellent chemical resistance for a variety of lab applications

Produced in 100K class clean room facility



Leak proof for ultimate protection and autoclavable*

Translucent visibility



USFDA approved virgin polypropylene of USP class VI

Available in PP, HDPE and LDPE



*All caps are made up of PP material

For more details visit us at www.accumaximum.com

AN ISO 9001 : 2015 & ISO 13485 : 2016 CERTIFIED MFG. COMPANY

*We also execute OEM projects. Specifications can be changed without notice for quality improvement.