

# ACCUMAX





VERSATILE, RELIABLE  
AND CONSISTENT

TECAN<sup>®</sup> TIPS

Volume: 10, 200 and 1000 ul

## TECAN® TIPS

The Accumax Robotic Tips for Tecan® have been carefully designed, produced, and validated to ensure reliable performance when used with LiHa Arms. These tips boast superior straightness and a low coefficient of variation, ensuring accurate and precise performance. They are compatible with the Freedom EVO® 100, 150, or 200 with Liquid Handling Arm, as well as the Freedom EVO® 75 8+1. Available in a variety of tip volumes, with or without filters, the Accumax Robotic Tips for Tecan® are the perfect choice for your automatic liquid handling needs.



## PRODUCT FEATURES



Constructed of premium-grade virgin polypropylene. Available in carbon-embedded conductive tips and clear tips

---



USP class VI resins

---



Validated on the actual instrument for which the tip is created

---

Manufactured in ISO-certified 9001 and 13485 facility

---



All tips are certified Endotoxin, Heavy Metals, RNase, DNase and DNA free

---




Nonsterile, sterile and filtered sterile formats.

---




## PRODUCT FEATURES



Polycarbonate (PC) tray with color coding to indicate different volumes


Polyethylene filter to prevent aerosol contamination



## Available In



Blister packs



Eco flip box

## Ordering Information

### Conductive tips

| Product Code   | Description   | Tray color | Per rack | Per unit      | Total |
|----------------|---|------------|----------|---------------|-------|
| ATTB-10-FL     | Accumax Tecan type 10µL tips<br>Standard-Conductive, Flip box                   | Grey       | 2x96     | 5x18 per pack | 17280 |
| ATTB-10-F-FL   | Accumax Tecan type 10µL tips<br>Standard Filter-Conductive, Flip Box            | Grey       | 2x96     | 5x18 per pack | 17280 |
| ATTB-200-FL    | Accumax Tecan type 200µL tips<br>Standard- conductive, Flip Box                 | Blue       | 2x96     | 5x18 per pack | 17280 |
| ATTB-200-F-FL  | Accumax Tecan type 200µL tips<br>Standard Filter-conductive, Flip Box           | Blue       | 2x96     | 5x18 per pack | 17280 |
| ATTB-1000-FL   | Accumax Tecan type 1000µL tips<br>Standard-conductive, Flip box                 | Red        | 2x96     | 5x10 per pack | 9600  |
| ATTB-1000-F-FL | Accumax Tecan type 1000µL tips<br>Standard-Filter- Conductive, Flip box         | Red        | 2x96     | 5x10 per pack | 9600  |
| ATTB-10-BL     | Accumax Tecan type 10µL tips<br>Standard-Conductive, 1-Blister pack             | Grey       | 1x96     | 24 Pack       | 2304  |
| ATTB-10-F-BL   | Accumax Tecan type 10µL tips<br>Standard Filter-Conductive, 1-<br>Blister pack  | Grey       | 1x96     | 24 Pack       | 2304  |
| ATTB-10-S-BL   | Accumax Tecan type 10µL tips<br>Sterile-Conductive, 1-Blister Pack              | Grey       | 1x96     | 24 Pack       | 2304  |
| ATTB-10-S-F-BL | Accumax Tecan type 10µL tips<br>Sterile-Filter-Conductive, 1-<br>Blister pack   | Grey       | 1x96     | 24 Pack       | 2304  |
| ATTB-200-BL    | Accumax Tecan type 200µL tips<br>Standard- conductive, 1-Blister pack           | Blue       | 1x96     | 24 Pack       | 2304  |
| ATTB-200-F-BL  | Accumax Tecan type 200µL tips<br>Standard Filter-conductive, 1-<br>Blister Pack | Blue       | 1x96     | 24 Pack       | 2304  |

| Product Code     | Description   | Tray color | Per rack | Per unit | Total |
|------------------|---|------------|----------|----------|-------|
| ATTB-200-S-BL    | Accumax tecan type 200µL tips<br>Sterile-conductive, 1-Blister Pack               | Blue       | 1x96     | 24 Pack  | 2304  |
| ATTB-200-S-F-BL  | Accumax Tecan type 200µL tips<br>Sterile-Filter-Conductive, 1-<br>Blister Pack    | Blue       | 1x96     | 24 Pack  | 2304  |
| ATTB-1000-BL     | Accumax Tecan type 1000µL tips<br>standard-conductive, 1-Blister pack             | Red        | 1x96     | 24 Pack  | 2304  |
| ATTB-1000-F-BL   | Accumax Tecan type 1000µL tips<br>Standard-Filter- Conductive, 1-<br>Blister pack | Red        | 1x96     | 24 Pack  | 2304  |
| ATTB-1000-S-BL   | Accumax Tecan type 1000µL tips<br>Sterile-conductive, 1-Blister pack              | Red        | 1x96     | 24 Pack  | 2304  |
| ATTB-1000-S-F-BL | Accumax Tecan type 1000µL tips<br>Sterile-Filter-Conductive, 1-<br>Blister pack   | Red        | 1x96     | 24 Pack  | 2304  |
| ATTB-50-FL       | Accumax Tecan type 50µL tips<br>Standard-Conductive, Flip box                     | Green      | 2x96     | 10 Pack  | 9600  |
| ATTB-50-F-FL     | Accumax Tecan type 50µL tips<br>Standard Filter-Conductive,<br>Flip box           | Green      | 2x96     | 10 Pack  | 9600  |
| ATTB-50-BL       | Accumax Tecan type 50µL tips<br>Standard-Conductive, Blister Pack                 | Green      | 1x96     | 24 Pack  | 2304  |
| ATTB-50-F-BL     | Accumax Tecan type 50µL tips<br>Standard Filter-Conductive,<br>Blister Pack       | Green      | 1x96     | 24 Pack  | 2304  |
| ATTB-50-S-BL     | Accumax Tecan type 50µL tips<br>Sterile-Conductive Blister Pack                   | Green      | 1x96     | 24 Pack  | 2304  |
| ATTB-50-S-F-BL   | Accumax Tecan type 50µL tips<br>Sterile-Filter-Conductive,<br>Blister Pack        | Green      | 1x96     | 24 Pack  | 2304  |

## Clear tips

| Product Code   | Description   | Tray color | Per rack | Per unit      | Total |
|----------------|---|------------|----------|---------------|-------|
| ATTC-10-FL     | Accumax Tecan type 10µL tips Standard- clear, Flip Box              | Grey       | 2x96     | 5x18 per pack | 17280 |
| ATTC-10-F-FL   | Accumax Tecan type 10µL tips Standard Filter- clear, Flip Box       | Grey       | 2x96     | 5x18 per pack | 17280 |
| ATTC-200-FL    | Accumax Tecan type 200µL tips Standard- clear, Flip Box             | Blue       | 2x96     | 5x18 per pack | 17280 |
| ATTC-200-F-FL  | Accumax Tecan type 200µL tips Standard Filter-clear, Flip Box       | Blue       | 2x96     | 5x18 per pack | 17280 |
| ATTC-1000-FL   | Accumax Tecan type 1000µL tips Standard-clear, Flip Box             | Red        | 2x96     | 5x10 per pack | 9600  |
| ATTC-1000-F-FL | Accumax Tecan type 1000µL tips Standard-Filter-clear, Flip Box      | Red        | 2x96     | 5x10 per pack | 9600  |
| ATTC-10-BL     | Accumax Tecan type 10µL tips Standard- clear, 1-Blister Pack        | Grey       | 1x96     | 24 Pack       | 2304  |
| ATTC-10-F-BL   | Accumax Tecan type 10µL tips Standard Filter- clear,1-Blister Pack  | Grey       | 1x96     | 24 Pack       | 2304  |
| ATTC-10-S-BL   | Accumax Tecan type 10µL tips Sterile-clear, 1-Blister Pack          | Grey       | 1x96     | 24 Pack       | 2304  |
| ATTC-10-S-F-BL | Accumax Tecan type 10µL tips Sterile-Filter-clear, 1-Blister pack   | Grey       | 1x96     | 24 Pack       | 2304  |
| ATTC-200-BL    | Accumax Tecan type 200µL tips Standard- clear,1-Blister pack        | Blue       | 1x96     | 24 Pack       | 2304  |
| ATTC-200-F-BL  | Accumax Tecan type 200µL tips Standard Filter-clear, 1-Blister pack | Blue       | 1x96     | 24 Pack       | 2304  |

| Product Code     | Description   | Tray color | Per rack | Per unit | Total |
|------------------|---|------------|----------|----------|-------|
| ATTC-200-S-BL    | Accumax Tecan type 200µL tips<br>Sterile-clear, 1-Blister pack          | Blue       | 1x96     | 24 Pack  | 2304  |
| ATTC-200-S-F-BL  | Accumax Tecan type 200µL tips<br>Sterile-Filter-clear, 1-Blister Pack   | Blue       | 1x96     | 24 Pack  | 2304  |
| ATTC-1000-BL     | Accumax Tecan type 1000µL tips<br>Standard-clear, 1-Blister Pack        | Red        | 1x96     | 24 Pack  | 2304  |
| ATTC-1000-F-BL   | Accumax Tecan type 1000µL tips<br>Standard-Filter-clear, 1-Blister pack | Red        | 1x96     | 24 Pack  | 2304  |
| ATTC-1000-S-BL   | Accumax Tecan type 1000µL tips<br>Sterile-clear, 1-Blister pack         | Red        | 1x96     | 24 Pack  | 2304  |
| ATTC-1000-S-F-BL | Accumax Tecan type 1000µL tips<br>Sterile-Filter-clear, 1-Blister pack  | Red        | 1x96     | 24 Pack  | 2304  |
| ATTC-10-FL       | Accumax Tecan type 50µL tips<br>Standard- clear Flip box                | Green      | 2x96     | 10 Pack  | 9600  |
| ATTC-50-F-FL     | Accumax Tecan type 50µL tips<br>Standard Filter- clear Flip box         | Green      | 2x96     | 10 Pack  | 9600  |
| ATTC-50-BL       | Accumax Tecan type 50µL tips<br>Standard- clear Flip box                | Green      | 1x96     | 24 Pack  | 2304  |
| ATTC-50-F-BL     | Accumax Tecan type 50µL tips<br>Standard- Filter clear Flip box         | Green      | 1x96     | 24 Pack  | 2304  |
| ATTC-50-S-BL     | Accumax Tecan type 50µL tips<br>Sterile – clear Flip box                | Green      | 1x96     | 24 Pack  | 2304  |
| ATTC-50-S-F-BL   | Accumax Tecan type 50µL tips<br>Sterile-Filter-clear Flip box           | Green      | 1x96     | 24 Pack  | 2304  |



# ACCUMAX

**Accumax Lab Devices Pvt. Ltd.**  
[www.accumaxlab.com](http://www.accumaxlab.com)

ISO 9001:2015  
ISO 13485:2016

\*We also execute OEM projects. Specifications can be changed without notice for quality improvement.