# **NEUATION**



# INNOVATION THAT MEETS YOUR REQUIREMENT

iSwix Jr.

Mini Vortex Mixer

#### **NEUATION**

## iSwix Jr.

Mini Vortex Mixer

When it comes to the usage of mixers in laboratories, what matters the most is their performance and reliability which could justify your investment. But, the issues remains that manufacturers are not aware of the struggles faced by users in their workplace and as such design machines that are frustrating to use and don't provide the kind of reliable performance that the user requires.

Here at Neuation, our R&D department takes the pain to understand the user's needs and their feedback to help us in designing a product which can bridge that gap between manufacturers and users.

As, a result we developed the iSwix Jr. which combines exceptional performance and rugged durability in a compact package. The vortex mixer is designed to offer thorough and quick mixing in a small body.



#### **NEUATION**

#### **Product Features**



Maximum load of 100 gm (3500 RPM - iSwix Jr.)



iSwix Jr. VT allows the user to set different speeds from 2500-3500 RPM



Orbital diameter of 4.5 mm



Brushless DC motor for ultra quiet performance and a maintenance-free long life



Inbuilt counter balance ensures stable operation for complete peace-of-mind Robust construction utilizes Die Cast Aluminium base and suction cups



Small footprint saves valuable bench space and can also be used inside hoods and cold rooms



IP43 compliant design makes it suitable for use



Safe Low Voltage Device



#### **Specifications**

iSwix Jr.		
Max RPM	3500 RPM	
Orbital diameter	4.5 mm	
Motor type	DC	
Variable speed	No	
Max tube diameter	30 mm	
Max Load Capacity	100 gm	
Press activation	Yes	
Speed Setting	Fixed	
Protection class	IP 43	
Dimensions (W X D X H)	100 x 100 x 67 mm	
Weight	0.55 kg	

### **Ordering Information**

iSwix Jr.		
Product Name	Cat #	Details
iSwix Jr.	JV-1110 EU/US/UK/AU	Mini Vortex Mixer with fix speed of 3500 RPM

