

ACCUMAX

**FINESSE &
CLARITY
PERSONIFIED**

ACCUMAX

PETRI

PLATES

ACCUMAX PETRI PLATES

Accumax produces high-quality Petri dish which is manufactured under strict aseptic conditions with USP class VI complaint polystyrene.

Designed specifically for microbiology applications. With a 90 mm diameter, they are perfect for use in automated plate filing systems. Vented and non-vented variants are available.

Vented dishes feature an elevated lid, promoting excellent gas exchange and evaporation. This design allows for optimal conditions for growth and incubation. On the other hand, non-vented dishes have a flat lid that prevents evaporation, encouraging a longer incubation period and facilitating anaerobic applications.



Product Features



Manufactured with USFDA
21 CFR compliant polystyrene



One compartment
with both vented and
non vented options



Fully automated
manufacturing and packaging



Endotoxin/Pyrogen free

Petri plates with
90 mm diameter



Produced in
Class 10K clean rooms



Raw materials are compliant
with ISO 1133 & 306 guidelines



Sterility Assurance level of
Sterile plates is 10^{-6} and that
of Aseptic plates is 10^{-3} SAL



Product Features



Optically clear petri plate



Two Packaging options:
10 plates/ bag & 20 plates/ bag



The plates have bevelled top and ridged bottom making it stackable and slipable for automated carousels.

Petri plates with excellent mechanical and heat resistance properties



Lot number, Sterility & expiration traceability in bag packaging



These plates are fully compatible with medium throughput carousel based media preparation systems.



Ordering Information

Petri Plates					
Cat No.	Product Description	Packing Type	Quantity / bag	Quantity bag / Case	Total Quantity per Case
AP-90V-10	Petriplates 90 mm diameter, Vented, Aseptic, bevelled top, raised ridge bottom, semi stackable and slipable	Bulk	10	50	500
AP-90NV-10	Petriplates 90 mm diameter, Non Vented, Aseptic, bevelled top, raised ridge bottom, semi stackable and slipable	Bulk	10	50	500
AP-90V-20	Petriplates 90 mm diameter, Vented, Aseptic, bevelled top, raised ridge bottom, semi stackable and slipable	Bulk	20	25	500
AP-90NV-20	Petriplates 90 mm diameter, Non Vented, Aseptic, bevelled top, raised ridge bottom, semi stackable and slipable	Bulk	20	25	500
AP-90V-10-S	Petriplates 90 mm diameter, Stack Size 10, Vented, Sterile, bevelled top, raised ridge bottom, semi stackable and slipable	Bulk	10	50	500
AP-90NV-10-S	Petriplates 90 mm diameter, Stack Size 10, Non Vented, Sterile, bevelled top, raised ridge bottom, semi stackable and slipable	Bulk	10	50	500
AP-90V-20-S	Petriplates 90 mm diameter, Stack Size 20, Vented, Sterile, bevelled top, raised ridge bottom, semi stackable and slipable	Bulk	20	25	500
AP-90NV-20-S	Petriplates 90 mm diameter, Stack Size 20, Non Vented, Sterile, bevelled top, raised ridge bottom, semi stackable and slipable	Bulk	20	25	500

ACCUMAX

Accumax Lab Devices Pvt. Ltd.
www.accumaxlab.com



ISO 9001:2015
ISO 13485:2016

*We also execute OEM projects. Specifications can be changed without notice for quality improvement.