

ACCUMAX PETRI PLATES



Manufactured with USFDA
21 CFR compliant polystyrene

Petri plates with
90 mm diameter



One compartment
with both vented and
non vented options

These plates are fully
compatible with medium
throughput carousel based
media preparation systems



Fully automated
manufacturing and packaging

Sterility Assurance level of
Sterile plates is 10^{-6} and that
of Aseptic plates is 10^{-3} SAL



Endotoxin/Pyrogen free

These Raised Ridge plates are
designed for easy stacking by hand.
They stack high and stay stacked
without slipping, sliding, or
scratching



Fully stackable with raised ridge
and optimally clear petri plate

Two Packaging options:
10 plates/ bag & 20 plates/
bag



For more details visit us at www.accumaximum.com