


# SOLUTION TO TIRED HANDS, BECAUSE PIPETTING SHOULDN'T HURT!

# ACCUMAX


## Accumax FAB LF Pipette



Effortless Low Force  
Tip Ejection



The Accumax Low Force  
Pipette forms a light, perfect  
seal between the pipette and  
tips, making tips virtually  
effortless to eject.



Designed with a cylindrical  
shaft and tip design that  
provides a reliable, consistent  
seal without excessive force.



Total hand comfort with  
lower spring forces.



Magnet  
assisted piston for  
consistently precise  
results



Innovative spring and  
seal design



Ergonomic design  
with proper weight  
distribution



Made using highly  
durable material for  
protection against  
chemical corrosion &  
physical shocks



UV resistant & Easy in-  
house calibration &  
tool free cleaning



Highly durable  
PVDF tip cone



Calibrated as per  
ISO 8655-6 In ISO 17025  
accredited lab



Fully  
autoclavable



For more details visit us at [www.accumaximum.com](http://www.accumaximum.com)

AN ISO 9001 : 2015 & ISO 13485 : 2016 CERTIFIED MFG. COMPANY

\*We also execute OEM projects. Specifications can be changed without notice for quality improvement.