

**A PERFECT
SOLUTION FOR
DIAGNOSTIC
SPECIMEN NEEDS**

**PETG
DIAGNOSTIC
VIAL**

MATERIAL - PETG, VOLUME RANGE : 5/10/20 ML



PETG DIAGNOSTIC VIAL

Material - PETG, Volume range : 5/10/20 ml

Accumax Polyethylene Terephthalate Glycol (PETG) Diagnostic vials with HDPE closures provide safe storage, transfer, and transport for a variety of important biological materials, including antibodies, proteins, nucleic acids, enzymes, substrates, reagents and buffers. They are ideal for clinical diagnostic laboratories.

Accumax PETG Diagnostic vials - the ultimate solution that meets all your biological and pharmaceutical regulatory standards. Crafted with precision, these vials are USP Class VI compliant, ensuring non-cytotoxic and non-pyrogenic properties for abnormal toxicity. Not only are they suitable for air shipment, but each manufacturing run undergoes rigorous testing to guarantee leak-proof performance.

What sets our PETG vials apart is their exclusive use of USP Class VI resin, specifically designed to meet the critical requirements of in-vitro fertilization media and all cell culture fluids.



Product Features



Manufactured in ISO approved clean room facility using USP Class VI Raw material



HDPE Cap closure with Polyethylene cap liner



Tested for low leachable and extractable



Tested to ensure leak proof performance

Excellent Transparency similar to glass



Gamma Sterile SAL10-6 product ensures maximum sterility



Non-Pyrogenic and Non-Cytotoxic



Ordering Information

PETG DIAGNOSTIC VIAL			
CAT Code	Description	Qty/Tray	Tray/Case
ADV-05-PG-S	5ml Diagnostic Vial with Cap, PETG, Sterile,Tray	20	5 Tray/100
ADV-10-PG-S	10ml Diagnostic Vial with Cap, PETG, Sterile, Tray	20	5 Tray/100
ADV-20-PG-S	20ml Diagnostic Vial with Cap, PETG, Sterile, Tray	20	5 Tray/100
ADV-05-PG-S-B	5ml Diagnostic Vial with Cap, PETG, Tray, Sterile, Bulk	100	5 Tray/500
ADV-10-PG-S-B	10ml Diagnostic Vial with Cap, PETG, Tray, Sterile, Bulk	100	5 Tray/500
ADV-20-PG-S-B	20ml Diagnostic Vial with Cap, PETG, Tray, Sterile, Bulk	100	5 Tray/500

ACCUMAX

Accumax Lab Devices Pvt. Ltd.
www.accumaximum.com

ISO 9001:2015
ISO 13485:2016

*We also execute OEM projects. Specifications can be changed without notice for quality improvement.