



Opt for GreenMAX and contribute
to sustainable environment


Accumax GreenMAX Tip Racks

Minimum Plastic – Maximum Environment!
Available for 10ul, 200ul, 1000ul and
1250ul Tips.

Accumax **GreenMAX** Tip Racks are the way to go if you're looking for innovative, sustainable lab consumables. Not only are they eco-friendly, but they also come in robust cardboard racks. They prove to be the perfect solution for those looking for a sustainable alternative to plastic racks. These racks work seamlessly with both single-channel and multi-channel pipettes, making loading and unloading tips a breeze. From an environmental standpoint, Accumax **GreenMAX** Tip Racks goes significantly beyond any other Tip rack to reduce its carbon footprint.

GreenMAX Tip Racks

Sustainable Packaging



Reduces your plastic consumption by 60%




Guilt-free disposal



Works seamlessly with single and multi-channel pipettes

Advantages over Plastic racks



- ✓ Saves space with a compact format that can be easily crushed and disposed.
- ✓ Light weight material that facilitates lower cost of disposal and transport.
- ✓ By using 100% biodegradable paper material in our tip racks, we are able to reduce plastic consumption by 60% compared to traditional tip racks.
- ✓ Material used is built to last, providing durability that rivals plastic racks. You can trust that your tips will remain intact until the very last one is used.

- ✓ Our technology minimizes the risk of particles collecting on the surface of the tips, reducing the possibility of contamination and error.
- ✓ **GreenMAX** Tip racks facilitates easy and guilt free disposal after use.
- ✓ By using these racks, you can eliminate plastic burden associated with traditional pipette tip storage and operation.
- ✓ These racks are non-autoclavable and pre-sterilized for immediate use.
- ✓ Tips and Tray Quality remains unchanged



Join us in our efforts to reduce plastic waste and switch to **ACCUMAX GreenMAX** Tip Racks today



Ordering Information

GreenMax ACCUMAX UNIVERSAL TIPS

Graduated, Free from DNase/ RNase/ Pyrogen/ Endotoxin

Catalog Number	Description	Packing Type	Quantity/Case
GT-10-R	10µl Universal Grad tip	Rack	96 Tips X 10 Racks X 10Packs
GT-10-S-R	10µl Universal Grad tip, Sterile	Rack	96 Tips X 10 Racks X 10Packs
GT-10-S-L-R	10µl Universal Grad tip, Sterile, Low retention	Rack	96 Tips X 10 Racks X 10Packs
GT-10-F-R(AT-10-F-K)	10µl Universal Grad tip, Filter	Rack	96 Tips X 10 Racks X 10Packs
GT-10-S-F-R(AT-10-S-F-K)	10µl Universal Grad tip, Sterile, Filter	Rack	96 Tips X 10 Racks X 10Packs
GT-10-S-L-F-R(AT-10-S-L-F-K)	10µl Universal Grad tip, Sterile, Low retention, Filter	Rack	96 Tips X 10 Racks X 10Packs
GT-10XL-R	10µl Long Universal Grad tip	Rack	96 Tips X 10 Racks X 10Packs
GT-10XL-S-R	10µl Long Universal Grad tip, 20µl capacity Sterile	Rack	96 Tips X 10 Racks X 10Packs
GT-10XL-S-L-R	10µl Long Universal Grad tip, 20µl capacity Sterile, Low retention	Rack	96 Tips X 10 Racks X 10Packs
GT-10XL-F-R(AT-10XL-F-K)	10µl Long Universal Grad tip, 20µl capacity Filter	Rack	96 Tips X 10 Racks X 10Packs
GT-10XL-S-F-R(AT-10XL-S-F-K)	10µl Long Universal Grad tip, 20µl capacity Sterile, Filter	Rack	96 Tips X 10 Racks X 10Packs
GT-10XL-S-L-F-R(AT-10XL-S-L-F-K)	10µl Long Universal Grad tip, 20µl capacity Sterile, Low retention, Filter	Rack	96 Tips X 10 Racks X 10Packs
GT-10XN-R	10µl Narrow Universal Grad tip	Rack	96 Tips X 10 Racks X 10Packs
GT-10XN-S-R	10µl Narrow Universal Grad tip, 20µl capacity Sterile	Rack	96 Tips X 10 Racks X 10Packs
GT-10XN-S-L-R	10µl Narrow Universal Grad tip, 20µl capacity Sterile, Low retention	Rack	96 Tips X 10 Racks X 10Packs
GT-10XN-F-R (AT-10XN-F-K)	10µl Narrow Universal Grad tip, 20µl capacity Filter	Rack	96 Tips X 10 Racks X 10Packs
GT-10XN-S-F-R(AT-10XN-S-F-K)	10µl Narrow Universal Grad tip, 20µl capacity Sterile, Filter	Rack	96 Tips X 10 Racks X 10Packs
GT-10XN-S-L-F-R(AT-10XN-S-L-F-K)	10µl Narrow Universal Grad tip, 20µl capacity Sterile, Low retention, Filter	Rack	96 Tips X 10 Racks X 10Packs
GT-20-F-R(AT-20-F-K)	20µl Universal Grad tip, Filter	Rack	96 Tips X 10 Racks X 10Packs
GT-20-S-F-R(AT-20-S-F-K)	20µl Universal Grad tip, Sterile, Filter	Rack	96 Tips X 10 Racks X 10Packs
GT-20-S-L-F-R(AT-20-S-L-F-K)	20µl Universal Grad tip, Sterile, Low retention, Filter	Rack	96 Tips X 10 Racks X 10Packs
GT-50-F-R (AT-50-F-K)	50µl Universal Grad tip, Filter	Rack	96 Tips X 10 Racks X 10Packs
GT-50-S-F-R(AT-50-S-F-K)	50µl Universal Grad tip, Sterile, Filter	Rack	96 Tips X 10 Racks X 10Packs
GT-50-S-L-F-R(AT-50-S-L-F-K)	50µl Universal Grad tip, Sterile, Low retention, Filter	Rack	96 Tips X 10 Racks X 10Packs
GT-100-F-R(AT-100-F-K)	100µl Universal Grad tip, Filter	Rack	96 Tips X 10 Racks X 10Packs
GT-100-S-F-R(AT-100-S-F-K)	100µl Universal Grad tip, Sterile, Filter	Rack	96 Tips X 10 Racks X 10Packs
GT-100-S-L-F-R(AT-100-S-L-F-K)	100µl Universal Grad tip, Sterile, Low retention, Filter	Rack	96 Tips X 10 Racks X 10Packs
GT-200-R	200µl Universal Grad tip	Rack	96 Tips X 10 Racks X 10Packs
GT-200-S-R	200µl Universal Grad tip, Sterile	Rack	96 Tips X 10 Racks X 10Packs
GT-200-S-L-R	200µl Universal Grad tip, Sterile, Low retention	Rack	96 Tips X 10 Racks X 10Packs
GT-200-F-R(AT-200-F-K)	200µl Universal Grad tip, Filter	Rack	96 Tips X 10 Racks X 10Packs
GT-200-S-F-R(AT-200-S-F-K)	200µl Universal Grad tip, Sterile, Filter	Rack	96 Tips X 10 Racks X 10Packs
GT-200-S-L-F-R(AT-200-S-L-F-K)	200µl Universal Grad tip, Sterile, Low retention, Filter	Rack	96 Tips X 10 Racks X 10Packs
GT-200XL-R	200µl Long Universal Grad tip, 300ul capacity	Rack	96 Tips X 10 Racks X 10Packs
GT-200XL-S-R	200µl Long Universal Grad tip, 300ul capacity, Sterile	Rack	96 Tips X 10 Racks X 10Packs
GT-200XL-S-L-R	200µl Long Universal Grad tip, 300ul capacity, Sterile, Low retention	Rack	96 Tips X 10 Racks X 10Packs
GT-200XL-F-R(AT-200XL-F-K)	200µl Long Universal Grad tip, 300ul capacity, Filter	Rack	96 Tips X 10 Racks X 10Packs
GT-200XL-S-F-R(AT-200XL-S-F-K)	200µl Long Universal Grad tip, 300ul capacity, Sterile, Filter	Rack	96 Tips X 10 Racks X 10Packs
GT-200XL-S-L-F-R(AT-200XL-S-L-F-K)	200µl Long Universal Grad tip, 300ul capacity, Sterile, Low retention, Filter	Rack	96 Tips X 10 Racks X 10Packs

Ordering Information

Catalog Number	Description	Packing Type	Quantity/Case
GT-1000-R	1000µl Universal Grad tip	Rack	96 Tips X 6 Racks X 10Packs
GT-1000-S-R	1000µl Universal Grad tip, Sterile	Rack	96 Tips X 6 Racks X 10Packs
GT-1000-S-L-R	1000µl Universal Grad tip, Sterile, Low retention	Rack	96 Tips X 6 Racks X 10Packs
GT-1000-F-R (AT-1000-F-K)	1000µl Universal Grad tip, Filter	Rack	96 Tips X 6 Racks X 10Packs
GT-1000-S-F-R(AT-1000-S-F-K)	1000µl Universal Grad tip, Sterile, Filter	Rack	96 Tips X 6 Racks X 10Packs
GT-1000-S-L-F-R(AT-1000-S-L-F-K)	1000µl Universal Grad tip, Sterile, Low retention, Filter	Rack	96 Tips X 6 Racks X 10Packs
GT-1000XL-R	1250µl Universal Grad tip	Rack	96 Tips X 6 Racks X 10Packs
GT-1000XL-S-R	1250µl Universal Grad tip, Sterile	Rack	96 Tips X 6 Racks X 10Packs
GT-1000XL-S-L-R	1250µl Universal Grad tip, Sterile, Low retention	Rack	96 Tips X 6 Racks X 10Packs
GT-1000XL-F-R(AT-1000XL-F-K)	1250µl Universal Grad tip, Filter	Rack	96 Tips X 6 Racks X 10Packs
GT-1000XL-S-F-R(AT-1000XL-S-F-K)	1250µl Universal Grad tip, Sterile, Filter	Rack	96 Tips X 6 Racks X 10Packs
GT-1000XL-S-L-F-R(AT-1000XL-S-L-F-K)	1250µl Universal Grad tip, Sterile, Low retention, Filter	Rack	96 Tips X 6 Racks X 10Packs

ACCUMAX

Accumax Lab Devices Pvt. Ltd.
www.accumaximum.com

ISO 9001:2015
ISO 13485:2016

*We also execute OEM projects. Specifications can be changed without notice for quality improvement.