

**A PERFECT
SOLUTION FOR
DIAGNOSTIC
SPECIMEN NEEDS**

**PETG
DIAGNOSTIC
VIAL**

MATERIAL - PETG, VOLUME RANGE : 5/10/20 ML



PETG DIAGNOSTIC VIAL

Material - PETG, Volume range : 5/10/20 ml

Accumax Polyethylene Terephthalate Glycol (PETG) Diagnostic vials with HDPE closures provide safe storage, transfer, and transport for a variety of important biological materials, including antibodies, proteins, nucleic acids, enzymes, substrates, reagents and buffers. They are ideal for clinical diagnostic laboratories.

Accumax PETG Diagnostic vials - the ultimate solution that meets all your biological and pharmaceutical regulatory standards. Crafted with precision, these vials are USP Class VI compliant, ensuring non-cytotoxic and non-pyrogenic properties for abnormal toxicity. Not only are they suitable for air shipment, but each manufacturing run undergoes rigorous testing to guarantee leak-proof performance.

What sets our PETG vials apart is their exclusive use of USP Class VI resin, specifically designed to meet the critical requirements of in-vitro fertilization media and all cell culture fluids.



Product Features



Manufactured in ISO approved clean room facility using USP Class VI Raw material



HDPE Cap* closure with Polyethylene cap liner



Tested for low leachable and extractable



Tested to ensure leak proof performance

Excellent Transparency similar to glass



Gamma Sterile SAL10-6 product ensures maximum sterility



Non-Pyrogenic and Non-Cytotoxic



Ordering Information

PETG DIAGNOSTIC VIAL			
CAT Code	Description	Qty/Tray	Tray/Case
ADV-05-PG-S	5ml Diagnostic Vial with Grey Cap, PETG, Sterile, Tray	20	5 Tray/100
ADV-10-PG-S	10ml Diagnostic Vial with Grey Cap, PETG, Sterile, Tray	20	5 Tray/100
ADV-20-PG-S	20ml Diagnostic Vial with Grey Cap, PETG, Sterile, Tray	20	5 Tray/100
ADV-05-PG-S-B	5ml Diagnostic Vial with Grey Cap, PETG, Tray, Sterile, Bulk	100	5 Tray/500
ADV-10-PG-S-B	10ml Diagnostic Vial with Grey Cap, PETG, Tray, Sterile, Bulk	100	5 Tray/500
ADV-20-PG-S-B	20ml Diagnostic Vial with Grey Cap, PETG, Tray, Sterile, Bulk	100	5 Tray/500

ACCUMAX

Accumax Lab Devices Pvt. Ltd.
www.accumaximum.com

ISO 9001:2015
ISO 13485:2016

*We also execute OEM projects. Specifications can be changed without notice for quality improvement.