

# PETG DIAGNOSTIC VIAL

# ACCUMAX

A perfect solution for diagnostic specimen needs.

Accumax PETG Diagnostic vials provide safe storage, transfer, and transport for a variety of important biological materials, including antibodies, proteins, nucleic acids, enzymes, and substrates.


**Vial Material** - PETG

**Cap Type** - Screw Cap

**Cap Closure Material** - HDPE

**Cap Liner** - Polyethylene


**Available Volume** - 5,10,20ml



Manufactured in ISO approved clean room facility using USP Class VI Raw material, with low leachable and extractable.



Tested to ensure leak proof performance.



Excellent Transparency similar to Glass.



Gamma Sterile product ensures maximum sterility.



Non-Pyrogenic and Non-Cytotoxic

For more details visit us at [www.accumaximum.com](http://www.accumaximum.com)

AN ISO 9001 : 2015 & ISO 13485 : 2016 CERTIFIED MFG. COMPANY

\*We also execute OEM projects. Specifications can be changed without notice for quality improvement.