

# WHEN PURITY IS PARAMOUNT

## Accumax Heavy Metal Free Tubes

# ACCUMAX



Pure, Additive-free  
Natural Color Caps  
For Easy Labelling



Tubes And Caps  
Made With Virgin,  
100% Medical-grade Resins.



Certified below one part  
per billion (<1ppb) trace  
of heavy metals




Large writing area



Heavy Metals Free centrifuge tubes are specifically manufactured to make sure that common trace metals that can interfere with experiments are below an ICP-MS detection level of one part per billion (<1ppb).

Metals included in testing are:  
Boron, Aluminum, Phosphorous,  
Arsenic, Selenium, Cadmium, Lead,  
Chromium, Manganese, Iron, Nickel,  
Copper, Zinc.



For more details visit us at [www.accumaximum.com](http://www.accumaximum.com)