

A NEW DIMENSION TO ACCURACY


ACCUMAX

Accumax Serological Pipette



Manufactured From Virgin Polystyrene In A Clean Room

Negative Graduation Marks For Extra Capacity



Certified Free From DNase, RNase, Cytotoxin Free & Non-pyrogenic

Adhering To 10^{-6} SAL (Sterility Assurance Level) Quality Standards



Calibrated For Accurate Dispensing Within +/- 2% Volume Range

Filter Plug Prevents Both Fluid And Liquid Vapors Minimizing Contamination Of Pipettor And Sample



Color Coding For Easy Identification Of Pipette Size



For more details visit us at www.accumaxlab.com

AN ISO 9001 : 2015 & ISO 13485 : 2016 CERTIFIED MFG. COMPANY

*We also execute OEM projects. Specifications can be changed without notice for quality improvement.