


ΔACCUMAX



PERFECTION REDEFINED

SCREW CAP TUBES

Volume : 0.5/1.5/2.0 ml

## Accumax SCREW CAP TUBES

Volume : 0.5/1.5/2.0 ml

Accumax Screw Cap Tubes provide a clean, safe environment for your critical samples, ensuring maximum protection and reliability. Designed with versatility in mind, these tubes range from 0.5ml to 2 ml and are highly customizable to meet diverse laboratory needs. Choose from options like graduations for precise measurements, ribbed grips for easy handling, and conical or skirted bases for specific applications.

Crafted with premium-quality material, these tubes ensure exceptional performance, whether you're using them for centrifugation, sample transport, or long-term storage. They are especially ideal for sensitive applications, including the secure handling of synthesized oligonucleotides, enzymes, buffers, and other high-value biological samples.

Built to withstand rigorous lab conditions, our Screw Cap Tubes are an indispensable tool for researchers and professionals who demand reliability and precision in their workflows. Enhance your lab's efficiency and confidence with these innovative, high-quality solutions.



0.5 ml Self  
Standing



1.5 ml Conical



1.5 ml Self  
Standing



1.5 ml Self  
Standing Amber




2.0 ml Conical



2.0 ml Self  
Standing

## Product Features



Accumax Screw Cap tubes are ideal for enzymes and Buffer storage.

Manufactured using USP Class VI medical grade Polypropylene (PP) conforming to US FDA



These tubes come in three Volumes 0.5ml, 1.5 ml & 2 ml and Self Standing/Conical formats.

They can withstand RCF upto 20,000g & -196 degree freezing & Thawing at Room temperature.



They are designed with an infused 'O' ring with better sealability and leakproof stability in comparison to traditional O ring design.

They are RNase free, DNase free, Endotoxin Free, Human DNA free, PCR Inhibitor free, ATP Free



Infused O ring, caps and tubes are all manufactured in Class 100K clean room with no secondary assembly operation involved making it contamination free production.

# Ordering Information

## Screw Cap Tubes With Separate Caps

0.5 ml Tubes				
Cat. No.	Description	Bag Quantity	Quantity / Pack	Quantity / Case
AST-005-T	0.5ml self-standing screw cap tube without cap	100	10x100	5 Pack
AST-005-T-S	0.5ml self standing screw cap tube Sterile without cap	100	10x100	5 Pack
AST-005-L-T-S	0.5ml self standing screw cap Low retention tube Sterile without cap	100	10x100	5 Pack
AST-005-T-C	0.5ml self standing screw cap tube with cap in separate bag	100	10x100	5 Pack
AST-005-T-S-C	0.5ml self standing screw cap tube Sterile with cap in separate bag	100	10x100	5 Pack
AST-005-L-T-S-C	0.5ml self standing screw cap Low retention tube Sterile with cap in separate bag	100	10x100	5 Pack
1.5 ml Tubes				
Cat. No.	Description	Bag Quantity	Quantity / Pack	Quantity / Case
ASCT-015-T	1.5ml Conical screw cap tube without cap	100	10x100	5 Pack
AST-015-T	1.5ml Self Standing screw cap tube without cap	100	10x100	5 Pack
ASCT-015-T-S	1.5ml Conical screw cap tube Sterile without cap	100	10x100	5 Pack
AST-015-T-S	1.5ml Self Standing screw cap tube Sterile without cap	100	10x100	5 Pack
ASCT-015-L-T-S	1.5ml Conical screw cap Low retention tube Sterile without cap	100	10x100	5 Pack
AST-015-L-T-S	1.5ml Self Standing screw cap Low retention tube Sterile without cap	100	10x100	5 Pack
ASCT-015-T-C	1.5ml Conical screw cap tube with cap in separate bag	100	10x100	5 Pack
AST-015-T-C	1.5ml Self Standing screw cap tube with cap in separate bag	100	10x100	5 Pack
ASCT-015-T-S-C	1.5ml Conical screw cap tube Sterile with cap in separate bag	100	10x100	5 Pack
AST-015-T-S-C	1.5ml Self Standing screw cap tube Sterile with cap in separate bag	100	10x100	5 Pack
ASCT-015-L-T-S-C	1.5ml Conical screw cap Low retention tube Sterile with cap in separate bag	100	10x100	5 Pack
AST-015-L-T-S-C	1.5ml Self Standing screw cap Low retention tube Sterile with cap in separate bag	100	10x100	5 Pack
2 ml Tubes				
Cat. No.	Description	Bag Quantity	Quantity / Pack	Quantity / Case
ASCT-020-T	2ml Conical screw cap tube without cap	100	10x100	5 Pack
AST-020-T	2 ml Self Standing screw cap tube without cap	100	10x100	5 Pack

ASCT-020-T-S	2ml Conical screw cap tube Sterile without cap	100	10x100	5 Pack
AST-020-T-S	2 ml Self Standing screw cap tube Sterile without cap	100	10x100	5 Pack
ASCT-020-L-T-S	2ml Conical screw cap Low retention tube Sterile without cap	100	10x100	5 Pack
AST-020-L-T-S	2 ml Self Standing screw cap Low retention tube Sterile without cap	100	10x100	5 Pack
ASCT-020-T-C	2ml Conical screw cap tube with cap in separate bag	100	10x100	5 Pack
AST-020-T-C	2 ml Self Standing screw cap tube with cap in separate bag	100	10x100	5 Pack
ASCT-020-T-S-C	2ml Conical screw cap tube Sterile with cap in separate bag	100	10x100	5 Pack
AST-020-T-S-C	2 ml Self Standing screw cap tube Sterile with cap in separate bag	100	10x100	5 Pack
ASCT-020-L-T-S-C	2ml Conical screw cap Low retention tube Sterile with cap in separate bag	100	10x100	5 Pack
AST-020-L-T-S-C	2 ml Self Standing screw cap Low retention tube Sterile with cap in separate bag	100	10x100	5 Pack

### Amber tubes 1.5ml

Cat. No.	Description	Bag Quantity	Quantity / Pack	Quantity / Case
AST-015-A-T	1.5ml Amber colour Self Standing screw cap tube without cap	100	10x100	5 Pack
AST-015-A-T-S	1.5ml Amber colour Self Standing screw cap tube Sterile without cap	100	10x100	5 Pack
AST-015-L-A-T-S	1.5ml Amber colour Self Standing screw cap Low retention tube Sterile without cap	100	10x100	5 Pack
AST-015-A-T-C	1.5ml Amber colour Self Standing screw cap tube with cap in separate bag	100	10x100	5 Pack
AST-015-A-T-S-C	1.5ml Amber colour Self Standing screw cap tube Sterile with cap in separate bag	100	10x100	5 Pack
AST-015-L-A-T-S-C	1.5ml Amber colour Self Standing screw cap Low retention tube Sterile with cap in separate bag	100	10x100	5 Pack

### Tube Caps

Cat. No.	Description	Bag Quantity	Quantity / Pack	Quantity / Case
AST-C-R	Accumax Universal screw caps Red colour	100	10x100	5 Pack
AST-C-Y	Accumax Universal screw caps Yellow colour	100	10x100	5 Pack
AST-C-G	Accumax Universal screw caps Green colour	100	10x100	5 Pack
AST-C-B	Accumax Universal screw caps Blue colour	100	10x100	5 Pack
AST-C-A	Accumax Universal screw caps Amber colour	100	10x100	5 Pack
AST-C-S-R	Accumax Universal screw caps Red colour, Sterile	100	10x100	5 Pack
AST-C-S-Y	Accumax Universal screw caps Yellow colour, Sterile	100	10x100	5 Pack
AST-C-S-G	Accumax Universal screw caps Green colour, Sterile	100	10x100	5 Pack
AST-C-S-B	Accumax Universal screw caps Blue colour, Sterile	100	10x100	5 Pack

AST-C-S-A	Accumax Universal screw caps Amber colour, Sterile	100	10x100	5 Pack
AST-C-Z	Accumax Universal screw caps Customize colour	100	10x100	5 Pack
AST-C-S-Z	Accumax Universal screw caps Customize colour, Sterile	100	10x100	5 Pack

ACCUMAX

Accumax Lab Devices Pvt. Ltd.  
www.accumaxlab.com

ISO 9001:2015  
ISO 13485:2016

\*We also execute OEM projects. Specifications can be changed without notice for quality improvement.