


# SCREW CAP TUBES

Volume : 0.5/1.5/2.0 ml

# ACCUMAX



Accumax Screw Cap tubes are ideal for enzymes and Buffer storage.



Manufactured using USP Class VI medical grade Polypropylene (PP) conforming to US FDA



These tubes come in three volumes 0.5ml, 1.5 ml & 2 ml and Self Standing/ conical formats.



They can withstand RCF upto 20,000g & -80 degree freezing & Thawing at Room temperature.



They are designed with an infused 'O' ring with better sealability and leakproof stability in comparison to traditional O ring design.



They are RNase free, DNase free, Endotoxin Free, Human DNA free, PCR Inhibitor free, ATP free



Infused O ring, caps and tubes are all manufactured in Class 100K clean room with no secondary assembly operation involved making it contamination free production.



0.5 ml Self Standing



1.5 ml Self Standing



1.5 ml Conical



1.5 ml Self Standing Amber



2.0 ml Conical



2.0 ml Self Standing

For more details visit us at [www.accumaxlab.com](http://www.accumaxlab.com)