

ACCUMAX







VERSATILE, RELIABLE
AND CONSISTENT

TECAN[®] TYPE TIPS

Volume: 10,50,200 and 1000 μ L

TECAN® TYPE TIPS

The Accumax Robotic Tips for Tecan® have been carefully designed, produced, and validated to ensure reliable performance when used with LiHa Arm® & Flexible Channel Arm (FCA)®. These tips boast superior straightness and a low coefficient of variation, ensuring accurate and precise performance. They are compatible with the Freedom EVO® 100, 150, or 200 with Liquid Handling Arm, the Freedom EVO® 75 8+1 as well as Fluent® FCA® (Flexible Channel Arm). Available in a variety of tip volumes, with or without filters, the Accumax Robotic Tips for Tecan® are the perfect choice for your automatic liquid handling needs.



PRODUCT FEATURES



Constructed of premium-grade virgin polypropylene. Available in carbon-embedded conductive tips and clear tips

Manufactured in ISO-certified 9001 and 13485 facility




USP class VI resins

All tips are certified Endotoxin, Heavy Metals, RNase, DNase and DNA free




Validated on the actual instrument for which the tip is created

Nonsterile, sterile and filtered sterile formats.




Polycarbonate (PC) tray with colorcoding to indicate volumes

Polyethylene filter to prevent aerosol contamination




Available In



BLISTER PACK



ANSI/SLAS RACK



ECO FLIP BOX

Ordering Information

Conductive tips

Product Code	Description	Tray color	Quantity / Pack	Quantity/ Case	Total Quantity
ATTB-10-FL	Accumax Tecan type 10µL tips Standard-Conductive, Flip box	Grey	2 x 96 tips x 5 Racks	10 pack	9600
ATTB-10-F-FL	Accumax Tecan type 10µL tips Standard Filter-Conductive, Flip Box	Grey	2 x 96 tips x 5 Racks	10 pack	9600
ATTB-50-F-FL	Accumax Tecan type 50µL tips Standard Filter-Conductive, Flip box	Green	2 x 96 tips x 5 Racks	10 Pack	9600
ATTB-50-FL	Accumax Tecan type 50µL tips Standard-Conductive, Flip box	Green	2 x 96 tips x 5 Racks	10 Pack	9600
ATTB-200-FL	Accumax Tecan type 200µL tips Standard- conductive, Flip Box	Blue	2 x 96 tips x 5 Racks	10 pack	9600
ATTB-200-F-FL	Accumax Tecan type 200µL tips Standard Filter-conductive, Flip Box	Blue	2 x 96 tips x 5 Racks	10 pack	9600
ATTB-1000-FL	Accumax Tecan type 1000µL tips Standard-conductive, Flip box	Yellow	2 x 96 tips x 5 Racks	10 pack	9600
ATTB-1000-F-FL	Accumax Tecan type 1000µL tips Standard-Filter- Conductive, Flip box	Yellow	2 x 96 tips x 5 Racks	10 pack	9600
ATTB-10-BL	Accumax Tecan type 10µL tips Standard-Conductive, 1-Blister pack	Grey	1x96 tips	24 Pack	2304
ATTB-10-F-BL	Accumax Tecan type 10µL tips Standard Filter-Conductive, 1-Blister pack	Grey	1x96 tips	24 Pack	2304
ATTB-10-S-BL	Accumax Tecan type 10µL tips Sterile-Conductive, 1-Blister Pack	Grey	1x96 tips	24 Pack	2304
ATTB-10-S-F-BL	Accumax Tecan type 10µL tips Sterile-Filter-Conductive, 1-Blister pack	Grey	1x96 tips	24 Pack	2304

Product Code	Description	Tray color	Quantity / Pack	Quantity/ Case	Total Quantity
ATTB-10-R-RF	Accumax Tecan type 10 µL Conductive Tips - ANSI/SLAS, Non-filtered, 1-Blister Pack	Grey	96	24	2304
ATTB-10-F-R-RF	Accumax Tecan type 10 µL Conductive Tips - ANSI/SLAS, Filtered, 1-Blister Pack	Grey	96	24	2304
ATTB-10-S-R	Accumax Tecan type 10 µL Conductive Tips - ANSI/SLAS, Non-filtered, Sterile, Single Stack	Grey	96	40	3840
ATTB-10-S-F-R	Accumax Tecan type 10 µL Conductive Tips - ANSI/SLAS, Filtered, Sterile, Single Stack	Grey	96	40	3840
ATTB-50-BL	Accumax Tecan type 50µL tips Standard-Conductive, 1-Blister Pack	Green	1x96 tips	24 Pack	2304
ATTB-50-F-BL	Accumax Tecan type 50µL tips Standard Filter-Conductive, 1-Blister pack	Green	1x96 tips	24 Pack	2304
ATTB-50-S-BL	Accumax Tecan type 50µL tips Sterile-Conductive, 1-Blister pack	Green	1x96 tips	24 Pack	2304
ATTB-50-S-F-BL	Accumax Tecan type 50µL tips Sterile-Filter-Conductive, 1-Blister pack	Green	1x96 tips	24 Pack	2304
ATTB-50-R- RF	Accumax Tecan type 50 µL Conductive Tips - ANSI/SLAS, Non-filtered, 1-Blister Pack	Green	96	24	2304
ATTB-50-F-R-RF	Accumax Tecan type 50 µL Conductive Tips - ANSI/SLAS, Filtered, 1-Blister Pack	Green	96	24	2304
ATTB-50-S-R	Accumax Tecan type 50 µL Conductive Tips - ANSI/SLAS, Non-filtered, Sterile, Single Stack	Green	96	40	3840
ATTB-50-S-F-R	Accumax Tecan type 50 µL Conductive Tips - ANSI/SLAS, Filtered, Sterile, Single Stack	Green	96	40	3840

Product Code	Description	Tray color	Quantity / Pack	Quantity/ Case	Total Quantity
ATTB-200-BL	Accumax Tecan type 200µL tips Standard- conductive, 1-Blister pack	Blue	1x96 tips	24 Pack	2304
ATTB-200-F-BL	Accumax Tecan type 200µL tips Standard Filter-conductive, 1-Blister Pack	Blue	1x96 tips	24 Pack	2304
ATTB-200-S-BL	Accumax tecan type 200µL tips Sterile-conductive, 1-Blister Pack	Blue	1x96 tips	24 Pack	2304
ATTB-200-S-F-BL	Accumax Tecan type 200µL tips Sterile-Filter-Conductive, 1-Blister Pack	Blue	1x96 tips	24 Pack	2304
ATTB-200-R-RF	Accumax Tecan type 200 µL Conductive Tips - ANSI/SLAS, Non-filtered, 1-Blister Pack	Blue	96	24	2304
ATTB-200-F-R-RF	Accumax Tecan type 200 µL Conductive Tips - ANSI/SLAS, Filtered, 1-Blister Pack	Blue	96	24	2304
ATTB-200-S-R	Accumax Tecan type 200 µL Conductive Tips - ANSI/SLAS, Non-filtered, Sterile, Single Stack	Blue	96	40	3840
ATTB-200-S-F-R	Accumax Tecan type 200 µL Conductive Tips - ANSI/SLAS, Filtered, Sterile, Single Stack	Blue	96	40	3840
ATTB-1000-BL	Accumax Tecan type 1000µL tips standard-conductive, 1-Blister pack	Yellow	1x96 tips	24 Pack	2304
ATTB-1000-F-BL	Accumax Tecan type 1000µL tips Standard-Filter- Conductive, 1-Blister pack	Yellow	1x96 tips	24 Pack	2304
ATTB-1000-S-BL	Accumax Tecan type 1000µL tips Sterile-conductive, 1-Blister pack	Yellow	1x96 tips	24 Pack	2304
ATTB-1000-S-F-BL	Accumax Tecan type 1000µL tips Sterile-Filter-Conductive, 1-Blister pack	Yellow	1x96 tips	24 Pack	2304

Product Code	Description	Tray color	Quantity / Pack	Quantity/ Case	Total Quantity
ATTB-1000-R-RF	Accumax Tecan type 1000 µL Conductive Tips - ANSI/SLAS, Non-filtered, 1-Blister Pack	Yellow	96	24	2304
ATTB-1000-F-R-RF	Accumax Tecan type 1000 µL Conductive Tips - ANSI/SLAS, Filtered, 1-Blister Pack	Yellow	96	24	2304
ATTB-1000-S-R	Accumax Tecan type 1000 µL Conductive Tips - ANSI/SLAS, Non-filtered, Sterile, Single Stack	Yellow	96	40	3840
ATTB-1000-S-F-R	Accumax Tecan type 1000 µL Conductive Tips - ANSI/SLAS, Filtered, Sterile, Single Stack	Yellow	96	40	3840

Clear tips

Product Code	Description	Tray color	Quantity / Pack	Quantity/ Case	Total Quantity
ATTC-10-FL	Accumax Tecan type 10µL tips Standard- clear, Flip Box	Grey	2 x 96 tips x 5 Racks	10 Pack	9600
ATTC-10-F-FL	Accumax Tecan type 10µL tips Standard Filter- clear, Flip Box	Grey	2 x 96 tips x 5 Racks	10 Pack	9600
ATTC-50-FL	Accumax Tecan type 50µL tips Standard- clear Flip box	Green	2 x 96 tips x 5 Racks	10 Pack	9600
ATTC-50-F-FL	Accumax Tecan type 50µL tips Standard Filter- clear Flip box	Green	2 x 96 tips x 5 Racks	10 Pack	9600
ATTC-200-FL	Accumax Tecan type 200µL tips Standard- clear, Flip Box	Blue	2 x 96 tips x 5 Racks	10 pack	9600
ATTC-200-F-FL	Accumax Tecan type 200µL tips Standard Filter-clear, Flip Box	Blue	2 x 96 tips x 5 Racks	10 pack	9600
ATTC-1000-FL	Accumax Tecan type 1000µL tips Standard-clear, Flip Box	Yellow	2 x 96 tips x 5 Racks	10 pack	9600

Product Code	Description	Tray color	Quantity / Pack	Quantity/ Case	Total Quantity
ATTC-1000-F-FL	Accumax Tecan type 1000µL tips Standard-Filter-clear, Flip Box	Yellow	2 x 96 tips x 5 Racks	10 pack	9600
ATTC-10-BL	Accumax Tecan type 10µL tips Standard- clear, 1-Blister Pack	Grey	1 x 96 tips	24 Pack	2304
ATTC-10-F-BL	Accumax Tecan type 10µL tips Standard Filter- clear,1-Blister Pack	Grey	1 x 96 tips	24 Pack	2304
ATTC-10-S-BL	Accumax Tecan type 10µL tips Sterile-clear, 1-Blister Pack	Grey	1 x 96 tips	24 Pack	2304
ATTC-10-S-F-BL	Accumax Tecan type 10µL tips Sterile-Filter-clear, 1-Blister pack	Grey	1 x 96 tips	24 Pack	2304
ATTC-10-R-RF	Accumax Tecan type 10 µL Clear Tips - ANSI/SLAS, Non-filtered, 1-Blister Pack	Grey	96	24	2304
ATTC-10-F-R-RF	Accumax Tecan type 10 µL Clear Tips - ANSI/SLAS,Filtered, 1-Blister Pack	Grey	96	24	2304
ATTC-10-S-R	Accumax Tecan type 10 µL Clear Tips - ANSI/SLAS, Non-filtered, Sterile, Single Stack	Grey	96	40	3840
ATTC-10-S-F-R	Accumax Tecan type 10 µL Clear Tips - ANSI/SLAS,Filtered, Sterile, Single Stack	Grey	96	40	3840
ATTC-50-BL	Accumax Tecan type 50µL tips Standard- clear,1-Blister pack	Green	1 x 96 tips	24 Pack	2304
ATTC-50-F-BL	Accumax Tecan type 50µL tips Standard- Filter clear,1-Blister pack	Green	1 x 96 tips	24 Pack	2304
ATTC-50-S-BL	Accumax Tecan type 50µL tips Sterile – clear,1-Blister pack	Green	1 x 96 tips	24 Pack	2304
ATTC-50-S-F-BL	Accumax Tecan type 50µL tips Sterile-Filter-clear,1-Blister pack	Green	1 x 96 tips	24 Pack	2304

Product Code	Description	Tray color	Quantity / Pack	Quantity/ Case	Total Quantity
ATTC-50-R- RF	Accumax Tecan type 50 µL ClearTips - ANSI/SLAS, Non-filtered, 1-Blister Pack	Green	96	24	2304
ATTC-50-F-R-RF	Accumax Tecan type 50 µL Clear Tips - ANSI/SLAS, Filtered, 1-Blister Pack	Green	96	24	2304
ATTC-50-S-R	Accumax Tecan type 50 µL Clear Tips - ANSI/SLAS,Non-filtered, Sterile, Single Stack	Green	96	40	3840
ATTC-50-S-F-R	Accumax Tecan type 50 µL Clear Tips - ANSI/SLAS, Filtered, Sterile, Single Stack	Green	96	40	3840
ATTC-200-BL	Accumax Tecan type 200µL tips Standard- clear,1-Blister pack	Blue	1 x 96 tips	24 Pack	2304
ATTC-200-F-BL	Accumax Tecan type 200µL tips Standard Filter-clear, 1-Blister pack	Blue	1 x 96 tips	24 Pack	2304
ATTC-200-S-BL	Accumax Tecan type 200µL tips Sterile-clear, 1-Blister pack	Blue	1 x 96 tips	24 Pack	2304
ATTC-200-S-F-BL	Accumax Tecan type 200µL tips Sterile-Filter-clear, 1-Blister Pack	Blue	1 x 96 tips	24 Pack	2304
ATTC-200-R-RF	Accumax Tecan type 200 µL Clear Tips - ANSI/SLAS, Non-filtered, 1-Blister Pack	Blue	96	24	2304
ATTC-200-F-R-RF	Accumax Tecan type 200 µL Clear Tips - ANSI/SLAS, Filtered, 1-Blister Pack	Blue	96	24	2304
ATTC-200-S-R	Accumax Tecan type 200 µL Clear Tips - ANSI/SLAS, Non-filtered, Sterile, Single Stack	Blue	96	40	3840
ATTC-200-S-F-R	Accumax Tecan type 200 µL Clear Tips - ANSI/SLAS, Filtered, Sterile, Single Stack	Blue	96	40	3840

Product Code	Description	Tray color	Quantity / Pack	Quantity/ Case	Total Quantity
ATTC-1000-BL	Accumax Tecan type 1000µL tips Standard-clear, 1-Blister Pack	Yellow	1 x 96 tips	24 Pack	2304
ATTC-1000-F-BL	Accumax Tecan type 1000µL tips Standard-Filter-clear, 1-Blister pack	Yellow	1 x 96 tips	24 Pack	2304
ATTC-1000-S-BL	Accumax Tecan type 1000µL tips Sterile-clear, 1-Blister pack	Yellow	1 x 96 tips	24 Pack	2304
ATTC-1000-S-F-BL	Accumax Tecan type 1000µL tips Sterile-Filter-clear, 1-Blister pack	Yellow	1 x 96 tips	24 Pack	2304
ATTC-1000-R-RF	Accumax Tecan type 1000 µL Clear Tips - ANSI/SLAS, Non-filtered, 1-Blister Pack	Yellow	96	24	2304
ATTC-1000-F-R-RF	Accumax Tecan type 1000 µL Clear Tips - ANSI/SLAS, Filtered, 1-Blister Pack	Yellow	96	24	2304
ATTC-1000-S-R	Accumax Tecan type 1000 µL Clear Tips - ANSI/SLAS, Non-filtered, Sterile, Single Stack	Yellow	96	40	3840
ATTC-1000-S-F-R	Accumax Tecan type 1000 µL Clear Tips - ANSI/SLAS, Filtered, Sterile, Single Stack	Yellow	96	40	3840

ACCUMAX

Accumax Lab Devices Pvt. Ltd.
www.accumaxlab.com

ISO 9001:2015
ISO 13485:2016

*We also execute OEM projects. Specifications can be changed without notice for quality improvement.